Package 'SeuratObject'

November 16, 2023

```
Type Package
Title Data Structures for Single Cell Data
Version 5.0.1
Description Defines S4 classes for single-cell genomic data and associated
     information, such as dimensionality reduction embeddings, nearest-neighbor
     graphs, and spatially-resolved coordinates. Provides data access methods and
     R-native hooks to ensure the Seurat object is familiar to other R users. See
     Satija R, Farrell J, Gennert D, et al (2015) <doi:10.1038/nbt.3192>,
     Macosko E, Basu A, Satija R, et al (2015) <doi:10.1016/j.cell.2015.05.002>,
     and Stuart T, Butler A, et al (2019) <doi:10.1016/j.cell.2019.05.031>,
     Hao Y, Hao S, et al (2021) <doi:10.1016/j.cell.2021.04.048> and
     Hao Y, et al (2023) <doi:10.1101/2022.02.24.481684> for
     more details.
URL https://mojaveazure.github.io/seurat-object/,
     https://github.com/mojaveazure/seurat-object
BugReports https://github.com/mojaveazure/seurat-object/issues
License MIT + file LICENSE
Encoding UTF-8
LazyData true
RoxygenNote 7.2.3
Additional_repositories https://bnprks.r-universe.dev
Depends R (>= 4.0.0),
     sp (>= 1.5.0)
Imports future,
     future.apply,
     grDevices,
     grid,
     Matrix (>= 1.6.3),
     methods,
     progressr,
     Rcpp (>= 1.0.5),
```

```
rlang (>= 0.4.7),
     stats,
     tools,
     utils,
     spam,
     lifecycle,
     generics
Suggests DelayedArray,
     fs (>= 1.5.2),
     ggplot2,
     HDF5Array,
     rmarkdown,
     testthat,
     BPCells,
     sf (>= 1.0.0)
Collate 'RcppExports.R'
     'zzz.R'
     'generics.R'
     'keymixin.R'
     'graph.R'
     'default.R'
     'assay.R'
     'logmap.R'
     'layers.R'
     'assay5.R'
     'centroids.R'
     'command.R'
     'compliance.R'
     'data.R'
     'jackstraw.R'
     'dimreduc.R'
     'segmentation.R'
     'molecules.R'
     'spatial.R'
     'fov.R'
     'neighbor.R'
     'seurat.R'
     'sparse.R'
     'utils.R'
LinkingTo Rcpp, RcppEigen
Enhances Seurat
Config/Needs/website pkgdown
R topics documented:
```

.DollarNames.SeuratCommand	7
AddMetaData	7
as.Centroids	8
as.Graph	9
as.list.SeuratCommand	10
as.matrix.LogMap	
as.Neighbor	11
as.Seurat	12
as.sparse	13
Assay-class	13
Assay-validity	14
Assay5-class	15
Assay5-validity	16
AssayData	16
Assays	
AttachDeps	
Boundaries	
CastAssay	
Cells	
CellsByIdentities	
CellsByImage	
Centroids-class	
Centroids-methods	
CheckGC	
CheckLayersName	
Command	
CreateAssay5Object	
CreateAssayObject	
CreateCentroids	
CreateDimReducObject	
CreateFOV	
CreateMolecules	
CreateSegmentation	
CreateSeuratObject	
Crop	
DefaultAssay	
DefaultDimReduc	
DefaultFOV	
DefaultLayer	
dim.Assay	
dim.Assay5	
dim.DimReduc	
dim.Seurat	
dimnames.Assay	
dimnames.Assay5	
dimnames.Seurat	
DimReduc-class	
DimReduc-validity	

Distances	49
droplevels.LogMap	50
Embeddings	51
FetchData	51
FilterObjects	53
FOV-class	
FOV-methods	
FOV-validity	
GetImage	
GetTissueCoordinates	
Graph-class	
HVFInfo	
Idents	
Images	
Index	
Indices	
intersect.LogMap	
IsGlobal	
IsMatrixEmpty	
IsNamedList	69
JackStrawData-class	70
JackStrawData-methods	70
JoinLayers	71
JS	
Key	
labels.LogMap	
LayerData	
Loadings	
LogMap	
LogMap-validity	
LogSeuratCommand	
merge.Assay	
merge.Assay5	
merge.DimReduc	
merge.Seurat	
Misc	
Molecules-class	
Molecules-methods	
names.Seurat	
Neighbor-class	
Neighbor-methods	90
Overlay	91
PackageCheck	92
pbmc_small	
print.DimReduc	
Project	
Radius	
RandomName	95

R	topics	documented:
---	--------	-------------

-5	
J	

RenameAssays	16
RenameCells	
RowMergeSparseMatrices	
SaveSeuratRds	
Segmentation-class	
Segmentation-methods	
set-if-null	
Seurat-class	
Seurat-validity)5
SeuratCommand-class)5
show,LogMap-method	16
Simplify	16
SpatialImage-class	17
SpatialImage-methods	17
split.Assay	
split.Assay5	
Stdev	
StitchMatrix	3
subset.Assay	4
subset.Assay5	
subset.DimReduc	5
subset.Seurat	6
Theta	7
Tool	
UpdateSeuratObject	
UpdateSlots	
Version	
WhichCells	1
[.Assay	
[.Assay5	
[.DimReduc	
[.SeuratCommand	
[[.Assay	
[[.Assay5	
[[.DimReduc	
[[.Seurat	
[[<-,Seurat	
[[<-,Seurat,NULL	
\$.Assay	
\$.Assay5	
\$.Seurat	
\$.SeuratCommand	5

SeuratObject-package SeuratObject: Data Structures for Single Cell Data

Description

Defines S4 classes for single-cell genomic data and associated information, such as dimensionality reduction embeddings, nearest-neighbor graphs, and spatially-resolved coordinates. Provides data access methods and R-native hooks to ensure the Seurat object is familiar to other R users. See Satija R, Farrell J, Gennert D, et al (2015) doi:10.1038/nbt.3192, Macosko E, Basu A, Satija R, et al (2015) doi:10.1016/j.cell.2015.05.002, and Stuart T, Butler A, et al (2019) doi:10.1016/j.cell.2019.05.031, Hao Y, Hao S, et al (2021) doi:10.1016/j.cell.2021.04.048 and Hao Y, et al (2023) doi:10.1101/2022.02.24.481684 for more details.

Author(s)

Maintainer: Paul Hoffman <seurat@nygenome.org> (ORCID)

Authors:

- Rahul Satija <rsatija@nygenome.org> (ORCID)
- Yuhan Hao <yhao@nygenome.org> (ORCID)
- Austin Hartman <ahartman@nygenome.org> (ORCID)
- Gesmira Molla <gmolla@nygenome.org> (ORCID)
- Andrew Butler <abutler@nygenome.org> (ORCID)
- Tim Stuart <tstuart@nygenome.org> (ORCID)

Other contributors:

- Madeline Kowalski <mkowalski@nygenome.org> (ORCID) [contributor]
- Saket Choudhary <schoudhary@nygenome.org> (ORCID) [contributor]
- Skylar Li <sli@nygenome.org> [contributor]
- Longda Jiang ljiang@nygenome.org> (ORCID) [contributor]
- Jeff Farrell < jfarrell@g.harvard.edu> [contributor]
- Shiwei Zheng <szheng@nygenome.org> (ORCID) [contributor]
- Christoph Hafemeister <chafemeister@nygenome.org> (ORCID) [contributor]

See Also

Useful links:

- https://mojaveazure.github.io/seurat-object/
- https://github.com/mojaveazure/seurat-object
- Report bugs at https://github.com/mojaveazure/seurat-object/issues

.DollarNames.SeuratCommand

Dollar-sign Autocompletion

Description

Autocompletion for \$ access on a SeuratCommand object

Usage

```
## S3 method for class 'SeuratCommand'
.DollarNames(x, pattern = "")
```

Arguments

x A SeuratCommand object

pattern A regular expression. Only matching names are returned.

Value

The parameter name matches for pattern

See Also

Command log object and interaction methods \$.SeuratCommand(), LogSeuratCommand(), SeuratCommand-class, [.SeuratCommand(), as.list.SeuratCommand()

AddMetaData

Add in metadata associated with either cells or features.

Description

Adds additional data to the object. Can be any piece of information associated with a cell (examples include read depth, alignment rate, experimental batch, or subpopulation identity) or feature (ENSG name, variance). To add cell level information, add to the Seurat object. If adding feature-level metadata, add to the Assay object (e.g. object[["RNA"]])

```
AddMetaData(object, metadata, col.name = NULL)
## S3 method for class 'Assay'
AddMetaData(object, metadata, col.name = NULL)
## S3 method for class 'Assay5'
```

as.Centroids

```
AddMetaData(object, metadata, col.name = NULL)
## S3 method for class 'Seurat'
AddMetaData(object, metadata, col.name = NULL)
```

Arguments

object An object

metadata A vector, list, or data.frame with metadata to add

col. name A name for meta data if not a named list or data.frame

Value

object with metadata added

Examples

```
cluster_letters <- LETTERS[Idents(object = pbmc_small)]
names(cluster_letters) <- colnames(x = pbmc_small)
pbmc_small <- AddMetaData(
  object = pbmc_small,
  metadata = cluster_letters,
  col.name = 'letter.idents'
)
head(x = pbmc_small[[]])</pre>
```

as.Centroids

Convert Segmentation Layers

Description

Convert Segmentation Layers

```
as.Centroids(x, nsides = NULL, radius = NULL, theta = NULL, ...)
as.Segmentation(x, ...)
## S3 method for class 'Segmentation'
as.Centroids(x, nsides = NULL, radius = NULL, theta = NULL, ...)
## S3 method for class 'Centroids'
as.Segmentation(x, ...)
```

as.Graph 9

Arguments

x An object
nsides The number of sides to represent cells/spots; pass Inf to plot as circles
radius Radius of shapes when plotting
theta Angle to adjust shapes when plotting
... Arguments passed to other methods

Value

```
as.Centroids: A Centroids object
as.Segmentation: A Segmentation object
```

as.Graph

Coerce to a Graph Object

Description

Convert a matrix (or Matrix) to a Graph object

Usage

```
as.Graph(x, ...)
## S3 method for class 'Matrix'
as.Graph(x, ...)
## S3 method for class 'matrix'
as.Graph(x, ...)
## S3 method for class 'Neighbor'
as.Graph(x, weighted = TRUE, ...)
```

Arguments

x The matrix to convert

... Ignored

weighted If TRUE, fill entries in Graph matrix with value from the nn.dist slot of the

Neighbor object

Value

A Graph object

See Also

Other graph: Graph-class

10 as.list.SeuratCommand

Examples

```
# converting sparse matrix
mat <- Matrix::rsparsematrix(nrow = 10, ncol = 10, density = 0.1)
rownames(x = mat) <- paste0("feature_", 1:10)
colnames(x = mat) <- paste0("cell_", 1:10)
g <- as.Graph(x = mat)

# converting dense matrix
mat <- matrix(data = 1:16, nrow = 4)
rownames(x = mat) <- paste0("feature_", 1:4)
colnames(x = mat) <- paste0("cell_", 1:4)
g <- as.Graph(x = mat)</pre>
```

as.list.SeuratCommand Coerce a SeuratCommand to a list

Description

Coerce a SeuratCommand to a list

Usage

```
## S3 method for class 'SeuratCommand'
as.list(x, complete = FALSE, ...)
```

Arguments

```
    x A SeuratCommand object
    complete Include slots besides just parameters (eg. call string, name, timestamp)
    Ignored
```

Value

A list with the parameters and, if complete = TRUE, the call string, name, and timestamp

See Also

Command log object and interaction methods \$.SeuratCommand(), .DollarNames.SeuratCommand(), LogSeuratCommand(), SeuratCommand-class, [.SeuratCommand()

Examples

```
cmd <- pbmc_small[["NormalizeData.RNA"]]
as.list(cmd)
as.list(cmd, complete = TRUE)</pre>
```

as.matrix.LogMap

as.matrix.LogMap

Coerce Logical Maps to Matrices

Description

Coerce a logical map to a matrix; this removes all logical map class capabilities from the object and returns a base-R matrix object

Usage

```
## S3 method for class 'LogMap'
as.matrix(x, ...)
```

Arguments

x A LogMap object
... Ignored

Value

A base-R matrix created from x

See Also

Logical map objects, validity, and interaction methods: LogMap-validity, LogMap, droplevels.LogMap(), intersect.LogMap(), labels.LogMap()

Examples

```
map <- LogMap(letters[1:10])
map[['obs']] <- c(1, 3, 7)
mat <- as.matrix(map)
mat
class(mat)</pre>
```

as.Neighbor

Coerce to a Neighbor Object

Description

Convert objects to Neighbor objects

12 as.Seurat

Usage

```
as.Neighbor(x, ...)
## S3 method for class 'Graph'
as.Neighbor(x, ...)
```

Arguments

x An object to convert to Neighbor

... Arguments passed to other methods

Value

A Neighbor object

as.Seurat

Coerce to a Seurat Object

Description

Convert objects to Seurat objects

Usage

```
as.Seurat(x, ...)
```

Arguments

x An object to convert to class Seurat

... Arguments passed to other methods

Value

A Seurat object generated from x

as.sparse 13

as.sparse

Cast to Sparse

Description

Convert dense objects to sparse representations

Usage

```
as.sparse(x, ...)
## S3 method for class 'data.frame'
as.sparse(x, row.names = NULL, ...)
## S3 method for class 'Matrix'
as.sparse(x, ...)
## S3 method for class 'matrix'
as.sparse(x, ...)
## S3 method for class 'ngCMatrix'
as.sparse(x, ...)
```

Arguments

x An object

. . . Arguments passed to other methods

row.names NULL or a character vector giving the row names for the data; missing values are

not allowed

Value

A sparse representation of the input data

Assay-class

The Assay Class

Description

The Assay object is the basic unit of Seurat; each Assay stores raw, normalized, and scaled data as well as cluster information, variable features, and any other assay-specific metadata. Assays should contain single cell expression data such as RNA-seq, protein, or imputed expression data.

14 Assay-validity

Slots

counts Unnormalized data such as raw counts or TPMs
data Normalized expression data
scale.data Scaled expression data
assay.orig Original assay that this assay is based off of. Used to track assay provenance
var.features Vector of features exhibiting high variance across single cells
meta.features Feature-level metadata

misc A named list of unstructured miscellaneous data

key A one-length character vector with the object's key; keys must be one or more alphanumeric characters followed by an underscore "_" (regex pattern "^[a-zA-Z][a-zA-Z0-9]*_\$")

See Also

```
v3 Assay object, validity, and interaction methods: $.Assay(), Assay-validity, CreateAssayObject(), [.Assay(), [[.Assay(), dim.Assay(), dimnames.Assay(), merge.Assay(), split.Assay(), subset.Assay()
```

Assay-validity

V3 Assay Validity

Description

Validation of Assay objects is handled by validObject

data Validation

blah

counts Validation

blah

scale.data Validation

blah

Feature-Level Meta Data Validation

blah

Variable Feature Validation

blah

Assay5-class 15

Key Validation

Keys must be a one-length character vector; a key must be composed of one of the following:

- An empty string (eg. "''") where nzchar() == 0
- An string composed of one or more alphanumeric values (both lower- and upper-case) that ends with an underscore ("_"); the first character must be a letter

Keys that are not empty strings are validated with the regex "^[a-zA-Z][a-zA-Z0-9]*_\$"

See Also

```
validObject
v3 Assay object, validity, and interaction methods: $.Assay(), Assay-class, CreateAssayObject(),
[.Assay(), [[.Assay(), dim.Assay(), dimnames.Assay(), merge.Assay(), split.Assay(),
subset.Assay()
```

Examples

```
rna <- pbmc_small[["RNA"]]
validObject(rna)</pre>
```

Assay5-class

The v5 Assay Object

Description

The v5 Assay is the typical Assay class used in **Seurat** v5; ...

Slots

- layers A named list containing expression matrices; each matrix should be a two-dimensional object containing some subset of cells and features defined in the cells and features slots. Cell and feature membership is recorded in the cells and features slots, respectively
- cells A logical mapping of cell names and layer membership; this map contains all the possible cells that this assay can contain. New layers must have some subset of cells present in this map
- features A <u>logical mapping</u> of feature names and layer membership; this map contains all the possible features that this assay can contain. New layers must have some subset of features present in this map
- default A one-length integer with the end index of the default layer; the default layer be all layers up to and including the layer at index default
- assay.orig Original assay that this assay is based off of; used to track assay provenance
- meta.data A data frame with feature-level meta data; should have the same number of rows as features
- misc A named list of unstructured miscellaneous data
- key A one-length character vector with the object's key; keys must be one or more alphanumeric characters followed by an underscore "_" (regex pattern "^[a-zA-Z][a-zA-Z0-9]*_\$")

16 AssayData

See Also

```
v5 Assay object, validity, and interaction methods: $.Assay5(), Assay5-validity, [.Assay5(), [[.Assay5(), dim.Assay5(), dimnames.Assay5(), merge.Assay5(), split.Assay5(), subset.Assay5()
```

Assay5-validity

V5 Assay Validity

Description

Validation of Assay5 objects is handled by validObject

Layer Validation

blah

Key Validation

Keys must be a one-length character vector; a key must be composed of one of the following:

- An empty string (eg. "''") where nzchar() == 0
- An string composed of one or more alphanumeric values (both lower- and upper-case) that ends with an underscore ("_"); the first character must be a letter

Keys that are not empty strings are validated with the regex "^[a-zA-Z][a-zA-Z0-9]*_\$"

See Also

```
validObject
```

```
v5 Assay object, validity, and interaction methods: $.Assay5(), Assay5-class, [.Assay5(), [[.Assay5(), dim.Assay5(), dimnames.Assay5(), merge.Assay5(), split.Assay5(), subset.Assay5()
```

AssayData

Get and Set Assay Data

Description

General accessor and setter functions for Assay objects. GetAssayData can be used to pull information from any of the expression matrices (eg. "counts", "data", or "scale.data"). SetAssayData can be used to replace one of these expression matrices

AssayData 17

Usage

```
GetAssayData(object, ...)
SetAssayData(object, layer, new.data, slot = deprecated(), ...)
## S3 method for class 'Seurat'
GetAssayData(object, assay = NULL, layer = NULL, slot = deprecated(), ...)
## S3 method for class 'Seurat'
SetAssayData(
  object,
 layer = "data",
 new.data,
  slot = deprecated(),
  assay = NULL,
)
## S3 method for class 'Assay'
GetAssayData(
  object,
  layer = c("data", "scale.data", "counts"),
  slot = deprecated(),
)
## S3 method for class 'Assay'
SetAssayData(
 object,
 layer = c("data", "scale.data", "counts"),
 new.data,
  slot = deprecated(),
)
```

Arguments

object	An object
	Arguments passed to other methods
layer	Name of layer to get or set
new.data	New assay data to add
slot	[Deprecated] Specific assay data to get or set
assay	Specific assay to get data from or set data for; defaults to the default assay

Value

GetAssayData: returns the specified assay data

18 Assays

SetAssayData: object with the assay data set

Lifecycle

[Superseded]

GetAssayData and SetAssayData have been superseded. To fetch expression matrices, use LayerData; to set expression data, use LayerData<-

Examples

```
# Get assay data from the default assay in a Seurat object
GetAssayData(object = pbmc_small, layer = "data")[1:5,1:5]
# Set an Assay layer through the Seurat object
count.data <- GetAssayData(object = pbmc_small[["RNA"]], layer = "counts")</pre>
count.data <- as.matrix(x = count.data + 1)</pre>
new.seurat.object <- SetAssayData(</pre>
    object = pbmc_small,
    layer = "counts",
    new.data = count.data,
    assay = "RNA"
)
# Get the data directly from an Assay object
GetAssayData(pbmc_small[["RNA"]], layer = "data")[1:5,1:5]
# Set an Assay layer directly
count.data <- GetAssayData(pbmc_small[["RNA"]], layer = "counts")</pre>
count.data \leftarrow as.matrix(x = count.data + 1)
new.assay <- SetAssayData(pbmc_small[["RNA"]], layer = "counts", new.data = count.data)</pre>
```

Assays

Query Specific Object Types

Description

List the names of Assay, DimReduc, Graph, Neighbor objects

```
Assays(object, ...)

Graphs(object, slot = NULL)

Neighbors(object, slot = NULL)

Reductions(object, slot = NULL)
```

AttachDeps 19

```
## S3 method for class 'Seurat'
Assays(object, slot = deprecated(), ...)
```

Arguments

object A Seurat object
... Ignored

slot Name of component object to return

Value

If slot is NULL, the names of all component objects in this Seurat object. Otherwise, the specific object specified

Examples

```
Assays(pbmc_small)
Graphs(pbmc_small)
Reductions(object = pbmc_small)
```

AttachDeps

Attach Required Packages

Description

Helper function to attach required packages. Detects if a package is already attached and if so, skips it. Should be called in .onAttach

Usage

```
AttachDeps(deps)
```

Arguments

deps

A character vector of packages to attach

Value

Invisibly returns NULL

Lifecycle

[Superseded]

AttachDeps has been superseded as of **SeuratObject** v5.0.0; as an alternative, list dependencies in the Depends section of DESCRIPTION

20 Boundaries

Examples

```
# Use in your .onAttach hook
if (FALSE) {
   .onAttach <- function(libname, pkgname) {
     AttachDeps(c("SeuratObject", "rlang"))
   }
}</pre>
```

Boundaries

Get, Set, and Query Segmentation Boundaries

Description

Get, Set, and Query Segmentation Boundaries

Usage

```
Boundaries(object, ...)

DefaultBoundary(object)

DefaultBoundary(object, ...) <- value

Molecules(object, ...)

## S3 method for class 'FOV'
Boundaries(object, ...)

## S3 method for class 'FOV'
DefaultBoundary(object)

## S3 replacement method for class 'FOV'
DefaultBoundary(object, ...) <- value

## S3 method for class 'FOV'
Molecules(object, ...)
```

Arguments

object An object

... Arguments passed to other methods

value The name of a segmentation boundary to set as default

CastAssay 21

Value

Boundaries: The names of all segmentation boundaries present within object

DefaultBoundary: The name of the default segmentation boundary

DefaultBoundary<-: object with the default segmentation boundary set to value

Molecules: The names of all molecule sets present within object

CastAssay

Cast Assay Layers

Description

Cast layers in v5 assays to other classes

Usage

```
CastAssay(object, to, ...)
## S3 method for class 'Assay5'
CastAssay(object, to, layers = NA, verbose = TRUE, ...)
```

Arguments

object	An object
to	Either a class name or a function that takes a layer and returns the same layer as a new class
• • •	If to is a function, arguments passed to to
layers	A vector of layers to cast; defaults to all layers
verbose	Show progress updates

Value

object with the layers cast to class specified by to

22 Cells

Cells

Cell and Feature Names

Description

Get the cell and feature names of an object

Usage

```
Cells(x, ...)
Features(x, ...)
## Default S3 method:
Cells(x, ...)
## S3 method for class 'Assay5'
Cells(x, layer = NULL, simplify = TRUE, ...)
## S3 method for class 'Assay5'
Features(x, layer = NULL, simplify = TRUE, ...)
## S3 method for class 'DimReduc'
Cells(x, ...)
## S3 method for class 'Neighbor'
Cells(x, ...)
```

Arguments

x An object

... Arguments passed to other methods

layer Layer to pull cells/features for; defaults to default layer; if NA, returns all cells

for the assay

simplify Simplify the cell/feature names into a single vector; if FALSE, separates each

cell/feature names by layer

Value

Cell: A vector of cell names

Features: A vector of feature names

See Also

```
dimnames.Assay5(), dimnames.Assay(), dimnames.Seurat()
```

CellsByIdentities 23

Examples

```
Cells(x = pbmc_small)
```

CellsByIdentities

Get cell names grouped by identity class

Description

Get cell names grouped by identity class

Usage

```
CellsByIdentities(object, idents = NULL, cells = NULL, return.null = FALSE)
```

Arguments

object A Seurat object

idents A vector of identity class levels to limit resulting list to; defaults to all identity

class levels

cells A vector of cells to grouping to

return.null If no cells are requested, return a NULL; by default, throws an error

Value

A named list where names are identity classes and values are vectors of cells belonging to that class

Examples

```
CellsByIdentities(object = pbmc_small)
```

 ${\tt CellsByImage}$

Get a vector of cell names associated with an image (or set of images)

Description

Get a vector of cell names associated with an image (or set of images)

```
CellsByImage(object, images = NULL, unlist = FALSE)
```

24 Centroids-class

Arguments

object Seurat object

images Vector of image names

unlist Return as a single vector of cell names as opposed to a list, named by image

name.

Value

A vector of cell names

Examples

```
## Not run:
CellsByImage(object = object, images = "slice1")
## End(Not run)
```

Centroids-class

The Centroids Class

Description

The Centroids Class

Slots

```
cells (character [n]) A vector of cell names; there should be as many cell names as there are points and no duplicate names
```

nsides (integer [1L]) The number of sides to draw when plotting centroids; must be either 0L for circles or greater than 3

radius (numeric [1L]) The radius of the shape when plotting the centroids

theta (numeric [1L]) The angle in degrees to adjust the shape when plotting the centroids

See Also

Centroids methods: Centroids-methods

Segmentation layer classes: Centroids-methods, Molecules-class, Molecules-methods, Segmentation-class, Segmentation-methods

Centroids-methods 25

Centroids-methods

Centroids Methods

Description

Methods for Centroids objects

```
## S3 method for class 'Centroids'
Cells(x, ...)
## S3 method for class 'Centroids'
GetTissueCoordinates(object, full = TRUE, ...)
## S3 method for class 'Centroids'
Radius(object)
## S3 method for class 'Centroids'
RenameCells(object, new.names = NULL, ...)
## S3 method for class 'Centroids'
Theta(object)
## S3 method for class 'Centroids'
is.finite(x)
## S3 method for class 'Centroids'
is.infinite(...)
## S3 method for class 'Centroids'
length(x)
## S3 method for class 'Centroids'
lengths(x, use.names = TRUE)
## S3 method for class 'Centroids'
subset(x, cells = NULL, ...)
## S4 method for signature 'Centroids, character, ANY, ANY'
x[i, j, ..., drop = TRUE]
## S4 method for signature 'Centroids, numeric, ANY, ANY'
x[i, j, ..., drop = TRUE]
## S4 method for signature 'Centroids'
show(object)
```

26 Centroids-methods

Arguments

x, object A Centroids object

... Arguments passed to other methods

full Expand the coordinates to the full polygon

new.names vector of new cell names

use.names Ignored

i, cells A vector of cells to keep; if NULL, defaults to all cells

j, drop Ignored

Details

GetTissueCoordinates: Get cell spatial coordinates

Radius: Get the centroid radius RenameCells: Update cell names

Theta: Get the offset angle

is.finite, is.infinite: Test to see if the centroids are circular or polygonal

length: Get the number of sides for the polygonal centroid lengths: Generate a run-length encoding of the cells present

subset, [: Subset a Centroids object to certain cells

show: Display an object summary to stdout

Value

GetTissueCoordinates: A data frame with three columns:

• "x": the x-coordinate

• "y": the y-coordinate

• "cell": the cell name

If full is TRUE, then each coordinate will indicate a vertex for the cell polygon; otherwise, each coordinate will indicate a centroid for the cell

Radius The radius of the centroids

RenameCells: object with the cells renamed to new.names

Theta: The offset angle in degrees

is.finite: TRUE if the centroids are polygonal, FALSE if circular

is.infinite: The opposite of is.finite

length: 0 if the centroids are circular, otherwise the number of sides of the polygonal centroid

lengths: An rle object for the cells

subset, [: x subsetted to the cells specified by cells/i

show: Invisibly returns NULL

CheckGC 27

See Also

Centroids-class

Segmentation layer classes: Centroids-class, Molecules-class, Molecules-methods, Segmentation-class, Segmentation-methods

CheckGC

Conditional Garbage Collection

Description

Call gc only when desired

Usage

```
CheckGC(option = "SeuratObject.memsafe")
```

Arguments

```
option ..
```

Value

Invisibly returns NULL

CheckLayersName

Check layers names for the input list

Description

Check layers names for the input list

Usage

```
CheckLayersName(matrix.list, layers.type = c("counts", "data"))
```

Arguments

```
matrix.list A list of matrices
```

layers.type layers type, such as counts or data

28 CreateAssay5Object

Command

Get SeuratCommands

Description

Pull information on previously run commands in the Seurat object.

Usage

```
Command(object, ...)
## S3 method for class 'Seurat'
Command(object, command = NULL, value = NULL, ...)
```

Arguments

object An object

... Arguments passed to other methods

command Name of the command to pull, pass NULL to get the names of all commands run

value Name of the parameter to pull the value for

Value

Either a SeuratCommand object or the requested parameter value

CreateAssay50bject

Create a v5 Assay object

Description

Create an Assay5 object from a feature expression matrix; the expected format of the matrix is features x cells

```
CreateAssay50bject(
  counts = NULL,
  data = NULL,
  min.cells = 0,
  min.features = 0,
  csum = NULL,
  fsum = NULL,
  ...
)
```

CreateAssayObject 29

Arguments

counts A two-dimensional expression matrix
data Optional prenormalized data matrix

min.cells Include features detected in at least this many cells; will subset the counts matrix

as well. To reintroduce excluded features, create a new object with a lower cutoff

fsum Function for calculating cell sums
Function for calculating feature sums
... Arguments passed to other methods

Value

An Assay5 object

Create Assay Object Create an Assay object

Description

Create an Assay object from a feature (e.g. gene) expression matrix. The expected format of the input matrix is features x cells.

Usage

```
CreateAssayObject(
  counts,
  data,
  min.cells = 0,
  min.features = 0,
  key = NULL,
  check.matrix = FALSE,
  ...
)
```

Arguments

counts Unnormalized data such as raw counts or TPMs
data Prenormalized data; if provided, do not pass counts

min.cells Include features detected in at least this many cells. Will subset the counts

matrix as well. To reintroduce excluded features, create a new object with a

lower cutoff

min. features Include cells where at least this many features are detected

key Optional key to initialize assay with

check.matrix Check counts matrix for NA, NaN, Inf, and non-integer values

... Arguments passed to as. sparse

30 CreateCentroids

Details

Non-unique cell or feature names are not allowed. Please make unique before calling this function.

Value

```
A Assay object
```

See Also

```
v3 Assay object, validity, and interaction methods: $.Assay(), Assay-class, Assay-validity, [.Assay(), [[.Assay(), dim.Assay(), dimnames.Assay(), merge.Assay(), split.Assay(), subset.Assay()
```

Examples

```
## Not run:
pbmc_raw <- read.table(
   file = system.file('extdata', 'pbmc_raw.txt', package = 'Seurat'),
   as.is = TRUE
)
pbmc_rna <- CreateAssayObject(counts = pbmc_raw)
pbmc_rna
## End(Not run)</pre>
```

CreateCentroids

Create a Centroids Objects

Description

Create a Centroids Objects

Usage

```
CreateCentroids(coords, nsides, radius, theta)
```

Arguments

coords	The coordinates of cell/spot centroids
nsides	The number of sides to represent cells/spots; pass Inf to plot as circles
radius	Radius of shapes when plotting
theta	Angle to adjust shapes when plotting

Value

A Centroids object

CreateDimReducObject Create a DimReduc object

Description

Create a DimReduc object

Usage

```
CreateDimReducObject(
  embeddings = new(Class = "matrix"),
  loadings = new(Class = "matrix"),
  projected = new(Class = "matrix"),
  assay = NULL,
  stdev = numeric(),
  key = NULL,
  global = FALSE,
  jackstraw = NULL,
  misc = list()
)
```

Arguments

embeddings	A matrix with the cell embeddings
loadings	A matrix with the feature loadings
projected	A matrix with the projected feature loadings
assay	Assay used to calculate this dimensional reduction
stdev	Standard deviation (if applicable) for the dimensional reduction
key	A character string to facilitate looking up features from a specific DimReduc
global	Specify this as a global reduction (useful for visualizations)
jackstraw	Results from the JackStraw function
misc	list for the user to store any additional information associated with the dimensional reduction

Value

A DimReduc object

See Also

```
Dimensional reduction object, validity, and interaction methods DimReduc-class, DimReduc-validity, [.DimReduc(), [[.DimReduc(), dim.DimReduc(), merge.DimReduc(), print.DimReduc(), subset.DimReduc()
```

32 CreateFOV

Examples

```
data <- GetAssayData(pbmc_small[["RNA"]], slot = "scale.data")
pcs <- prcomp(x = data)
pca.dr <- CreateDimReducObject(
  embeddings = pcs$rotation,
  loadings = pcs$x,
  stdev = pcs$sdev,
  key = "PC",
  assay = "RNA"
)</pre>
```

CreateF0V

Create Spatial Coordinates

Description

Create Spatial Coordinates

```
CreateFOV(coords, ...)
## S3 method for class 'Centroids'
CreateFOV(
  coords,
 molecules = NULL,
  assay = "Spatial",
 key = NULL,
 name = NULL,
  . . .
)
## S3 method for class 'data.frame'
CreateFOV(
  coords,
  type = c("segmentation", "centroids"),
  nsides = Inf,
  radius = NULL,
  theta = 0L,
 molecules = NULL,
  assay = "Spatial",
  key = NULL,
  name = NULL,
)
```

CreateMolecules 33

```
## S3 method for class 'list'
CreateFOV(coords, molecules = NULL, assay = "Spatial", key = NULL, ...)
## S3 method for class 'Segmentation'
CreateFOV(
   coords,
   molecules = NULL,
   assay = "Spatial",
   key = NULL,
   name = NULL,
   ...
)
```

Arguments

coords

type

... Arguments passed to other methods

Molecules

A data.frame with spatially-resolved molecule information or a Molecules object

Assay

Name of associated assay

Key

Key for these spatial coordinates

Name

When coords is a data.frame, Centroids, or Segmentation, name to store coordinates as

When providing a data.frame, specify if the coordinates represent a cell seg-

mentation or voxel centroids

nsides The number of sides to represent cells/spots; pass Inf to plot as circles

radius Radius of shapes when plotting
theta Angle to adjust shapes when plotting

Spatial coordinates

Value

A FOV object

See Also

FOV-class

Create Molecules Create a Molecules Object

Description

Create a Molecules Object

34 CreateSegmentation

Usage

```
CreateMolecules(coords, ...)
## S3 method for class 'data.frame'
CreateMolecules(coords, key = "", ...)
## S3 method for class 'Molecules'
CreateMolecules(coords, ...)
## S3 method for class 'NULL''
CreateMolecules(coords, ...)
```

Arguments

coords

Spatial coordinates for molecules; should be a data frame with three columns:

- "x": x-coordinates for each molecule
- "y": y-coordinates for each molecule
- "gene": gene name for each molecule

... Arguments passed to other methods

key A key to set for the molecules

Value

A Molecules object

CreateSegmentation

Create a Segmentation Objects

Description

Create a Segmentation Objects

Usage

```
CreateSegmentation(coords)
## S3 method for class 'data.frame'
CreateSegmentation(coords)
## S3 method for class 'Segmentation'
CreateSegmentation(coords)
```

Arguments

coords

The coordinates of cell segmentations

CreateSeuratObject 35

Value

A Segmentation object

CreateSeuratObject

Create a Seurat object

Description

Create a Seurat object from raw data

```
CreateSeuratObject(
  counts,
  assay = "RNA",
 names.field = 1,
 names.delim = "_",
 meta.data = NULL,
 project = "CreateSeuratObject",
)
## Default S3 method:
CreateSeuratObject(
  counts,
  assay = "RNA",
 names.field = 1L,
 names.delim = "_",
 meta.data = NULL,
 project = "SeuratProject",
 min.cells = 0,
 min.features = 0,
)
## S3 method for class 'Assay'
CreateSeuratObject(
  counts,
  assay = "RNA",
  names.field = 1L,
 names.delim = "_",
 meta.data = NULL,
 project = "SeuratProject",
)
```

36 CreateSeuratObject

```
## S3 method for class 'Assay5'
CreateSeuratObject(
  counts,
  assay = "RNA",
  names.field = 1L,
  names.delim = "_",
  meta.data = NULL,
  project = "SeuratProject",
  ...
)
```

Arguments

counts

features as rows or an Assay-derived object Name of the initial assay assay names.field For the initial identity class for each cell, choose this field from the cell's name. E.g. If your cells are named as BARCODE_CLUSTER_CELLTYPE in the input matrix, set names.field to 3 to set the initial identities to CELLTYPE. For the initial identity class for each cell, choose this delimiter from the cell's names.delim column name. E.g. If your cells are named as BARCODE-CLUSTER-CELLTYPE, set this to "-" to separate the cell name into its component parts for picking the relevant field. Additional cell-level metadata to add to the Seurat object. Should be a data.frame meta.data where the rows are cell names and the columns are additional metadata fields. Row names in the metadata need to match the column names of the counts matrix. project Project name for the Seurat object Arguments passed to other methods Include features detected in at least this many cells. Will subset the counts min.cells

Either a matrix-like object with unnormalized data with cells as columns and

matrix as well. To reintroduce excluded features, create a new object with a

Value

A Seurat object

min.features

Note

In previous versions (<3.0), this function also accepted a parameter to set the expression threshold for a 'detected' feature (gene). This functionality has been removed to simplify the initialization process/assumptions. If you would still like to impose this threshold for your particular dataset, simply filter the input expression matrix before calling this function.

Include cells where at least this many features are detected

Crop 37

Examples

```
## Not run:
pbmc_raw <- read.table(
   file = system.file('extdata', 'pbmc_raw.txt', package = 'Seurat'),
   as.is = TRUE
)
pbmc_small <- CreateSeuratObject(counts = pbmc_raw)
pbmc_small
## End(Not run)</pre>
```

Crop

Crop Coordinates

Description

Crop Coordinates

Usage

```
Crop(object, x = NULL, y = NULL, coords = c("plot", "tissue"), ...)
## S3 method for class 'FOV'
Crop(object, x = NULL, y = NULL, coords = c("plot", "tissue"), ...)
```

Arguments

object	An object
x, y	Range to crop x/y limits to; if NULL, uses full range of x/y
coords	Coordinate system to execute crop; choose from:
	• "plot": Coordinates as shown when plotting
	• "tissue": Coordinates from GetTissueCoordinates
	Arguments passed to other methods

Value

```
object cropped to the region specified by x and y
```

38 DefaultAssay

DefaultAssay

Default Assay

Description

Get and set the default assay

Usage

```
DefaultAssay(object, ...)
DefaultAssay(object, ...) <- value</pre>
## S3 method for class 'Graph'
DefaultAssay(object, ...)
## S3 replacement method for class 'Graph'
DefaultAssay(object, ...) <- value</pre>
## S3 method for class 'Assay'
DefaultAssay(object, ...)
## S3 replacement method for class 'Assay'
DefaultAssay(object, ...) <- value</pre>
## S3 method for class 'Assay5'
DefaultAssay(object, ...)
## S3 replacement method for class 'Assay5'
DefaultAssay(object, ...) <- value</pre>
## S3 method for class 'SeuratCommand'
DefaultAssay(object, ...)
## S3 method for class 'DimReduc'
DefaultAssay(object, ...)
## S3 replacement method for class 'DimReduc'
DefaultAssay(object, ...) <- value</pre>
## S3 method for class 'Seurat'
DefaultAssay(object, ...)
## S3 replacement method for class 'Seurat'
DefaultAssay(object, ...) <- value</pre>
```

DefaultDimReduc 39

Arguments

object An object

... Arguments passed to other methods value Name of assay to set as default

Value

 $\label{eq:default} \textbf{DefaultAssay: The name of the default assay}$

DefaultAssay<-: An object with the default assay updated

Examples

```
# Get current default assay
DefaultAssay(object = pbmc_small)

# Create dummy new assay to demo switching default assays
new.assay <- pbmc_small[["RNA"]]
Key(object = new.assay) <- "RNA2_"
pbmc_small[["RNA2"]] <- new.assay
# switch default assay to RNA2
DefaultAssay(object = pbmc_small) <- "RNA2"
DefaultAssay(object = pbmc_small)</pre>
```

DefaultDimReduc

Find the default DimReduc

Description

Searches for DimReducs matching "umap", "tsne", or "pca", case-insensitive, and in that order. Priority given to DimReducs matching the DefaultAssay or assay specified (eg. "pca" for the default assay weights higher than "umap" for a non-default assay)

Usage

```
DefaultDimReduc(object, assay = NULL)
```

Arguments

object A Seurat object

assay Name of assay to use; defaults to the default assay of the object

Value

The default DimReduc, if possible

DefaultFOV

Examples

```
DefaultDimReduc(pbmc_small)
```

DefaultF0V

Get and Set the Default FOV

Description

Get and Set the Default FOV

Usage

```
DefaultFOV(object, ...)
DefaultFOV(object, ...) <- value

## S3 method for class 'Seurat'
DefaultFOV(object, assay = NULL, ...)

## S3 replacement method for class 'Seurat'
DefaultFOV(object, assay = NA, ...) <- value</pre>
```

Arguments

object	A Seurat Object
	Arguments passed to other methods
value	The name of the FOV to set as the default
assay	Name of assay to get or set default FOV for; pass NA to get or set the global default ${\sf FOV}$

Value

```
DefaultFOV: The name of the default FOV

DefaultFOV<-: object with the default FOV set to value
```

DefaultLayer 41

DefaultLayer

Default Layer

Description

Get and set the default layer

Usage

```
DefaultLayer(object, ...)
DefaultLayer(object, ...) <- value

## S3 method for class 'Assay'
DefaultLayer(object, ...)

## S3 method for class 'Assay5'
DefaultLayer(object, ...)

## S3 replacement method for class 'Assay5'
DefaultLayer(object, ...) <- value</pre>
```

Arguments

object An object

... Arguments passed to other methods

value Name of layer to set as default

Value

DefaultLayer: The name of the default layer

DefaultLayer<-: An object with the default layer updated

dim.Assay

Feature and Cell Numbers

Description

Feature and Cell Numbers

Usage

```
## S3 method for class 'Assay'
dim(x)
```

42 dim.Assay5

Arguments

Х

An Assay object

Value

A two-length numeric vector with the total number of features and cells in x

See Also

```
v3 Assay object, validity, and interaction methods: $.Assay(), Assay-class, Assay-validity,
CreateAssayObject(), [.Assay(), [[.Assay(), dimnames.Assay(), merge.Assay(), split.Assay(),
subset.Assay()
```

Examples

```
rna <- pbmc_small[["RNA"]]
dim(rna)</pre>
```

dim.Assay5

Feature and Cell Numbers

Description

Feature and Cell Numbers

Usage

```
## S3 method for class 'Assay5'
dim(x)
```

Arguments

Х

An Assay5 object

Value

A two-length numeric vector with the total number of features and cells in x

See Also

```
v5 Assay object, validity, and interaction methods: $.Assay5(), Assay5-class, Assay5-validity, [.Assay5(), [[.Assay5(), dimnames.Assay5(), merge.Assay5(), split.Assay5(), subset.Assay5()
```

dim.DimReduc 43

dim.DimReduc

Dimensional Reduction Meta-Information

Description

Pull meta-information about cells and dimensions for a given dimensional reduction; cell meta-information is stored as row meta-information (eg. nrow, rownames) and dimension meta-information is stored as column meta-information (eg. ncol, colnames)

Usage

```
## S3 method for class 'DimReduc'
dim(x)

## S3 method for class 'DimReduc'
dimnames(x)

## S3 method for class 'DimReduc'
length(x)

## S3 method for class 'DimReduc'
names(x)
```

Arguments

Х

A DimReduc object

Value

```
dim: The number of cells (nrow) and dimensions (ncol) dimnames: The cell (row) and dimension (column) names length: The number of dimensions names: The dimension identifiers
```

See Also

Cells

```
Dimensional reduction object, validity, and interaction methods CreateDimReducObject(), DimReduc-class, DimReduc-validity, [.DimReduc(), [[.DimReduc(), merge.DimReduc(), print.DimReduc(), subset.DimReduc()
```

```
pca <- pbmc_small[["pca"]]
pca
dim(pca)</pre>
```

44 dim.Seurat

```
# nrow is number of cells
nrow(pca)

# rownames pulls cell names
head(rownames(pca))

# ncol and length are number of dimensions
ncol(pca)
length(pca)

# colnames and names pull dimension identifiers
head(colnames(pca))
head(names(pca))
```

dim.Seurat

Feature and Cell Numbers

Description

Feature and Cell Numbers

Usage

```
## S3 method for class 'Seurat'
dim(x)
```

Arguments

Χ

A Seurat object

Value

A two-length numeric vector with the total number of features and cells in x

See Also

```
Seurat object, validity, and interaction methods $.Seurat(), Seurat-class, Seurat-validity,
[[.Seurat(), [[<-,Seurat, NULL, [[<-,Seurat, dimnames.Seurat(), merge.Seurat(), names.Seurat(),
subset.Seurat()</pre>
```

```
# Get the number of features in an object
nrow(pbmc_small)
# Get the number of cells in an object
ncol(pbmc_small)
```

dimnames.Assay 45

dimnames.Assay

Assay-Level Feature and Cell Names

Description

Get and set feature and cell names in v5 Assays

Usage

```
## S3 method for class 'Assay'
dimnames(x)

## S3 replacement method for class 'Assay'
dimnames(x) <- value</pre>
```

Arguments

x An Assay object

value A two-length list where the first entry is the existing feature names for x and the

second entry is the *updated* cell names for x

Value

dimnames: A two-length list with the following values:

- A character vector will all features in x
- A character vector will all cells in x

dimnames<-: x with the cell names updated to those in value[[2L]]

See Also

```
v3 Assay object, validity, and interaction methods: $.Assay(), Assay-class, Assay-validity,
CreateAssayObject(), [.Assay(), [[.Assay(), dim.Assay(), merge.Assay(), split.Assay(),
subset.Assay()
Cells(), dimnames.Assay5(), dimnames.Seurat()
```

```
rna <- pbmc_small[["RNA"]]
# Feature and cell names can be acquired with `rownames` and `colnames`
head(rownames(rna))
head(colnames(rna))
# Cell names can be updated with `colnames<-`
colnames(rna)[1] <- "newcell"
head(colnames(rna))</pre>
```

46 dimnames.Seurat

dimnames.Assay5

Assay-Level Feature and Cell Names

Description

Get and set feature and cell names in v5 Assays

Usage

```
## S3 method for class 'Assay5'
dimnames(x)
## S3 replacement method for class 'Assay5'
dimnames(x) <- value</pre>
```

Arguments

x An Assay5 object

value A two-length list with updated feature and/or cells names

Value

dimnames: A two-length list with the following values:

- A character vector with all features in x
- A character vector with all cells in x

dimnames<-: x with the feature and/or cell names updated to value

See Also

```
v5 Assay object, validity, and interaction methods: $.Assay5(), Assay5-class, Assay5-validity, [.Assay5(), [[.Assay5(), dim.Assay5(), merge.Assay5(), split.Assay5(), subset.Assay5() Cells(), dimnames.Assay(), dimnames.Seurat()
```

dimnames.Seurat

Feature and Cell Names

Description

Get and set feature and cell inames in Seurat objects

DimReduc-class 47

Usage

```
## S3 method for class 'Seurat'
dimnames(x)

## S3 replacement method for class 'Seurat'
dimnames(x) <- value</pre>
```

Arguments

x A Seurat object

value A two-length list with updated feature and/or cells names

Value

dimnames: A two-length list with the following values:

- A character vector with all features in the default assay
- A character vector with all cells in x

dimnames<-: x with the feature and/or cell names updated to value

See Also

```
Seurat object, validity, and interaction methods $.Seurat(), Seurat-class, Seurat-validity,
[[.Seurat(), [[<-, Seurat, NULL, [[<-, Seurat, dim.Seurat(), merge.Seurat(), names.Seurat(),
subset.Seurat()
Cells(), dimnames.Assay5(), dimnames.Assay()</pre>
```

Examples

```
# Get the feature names of an object
head(rownames(pbmc_small))

# Get the cell names of an object
head(colnames(pbmc_small))

colnames(pbmc_small)[1] <- "newcell"
head(colnames(pbmc_small))</pre>
```

DimReduc-class

The Dimensional Reduction Class

Description

The DimReduc object stores a dimensionality reduction taken out in Seurat; each DimReduc consists of a cell embeddings matrix, a feature loadings matrix, and a projected feature loadings matrix.

48 DimReduc-validity

Slots

```
cell.embeddings Cell embeddings matrix (required)

feature.loadings Feature loadings matrix (optional)

feature.loadings.projected Projected feature loadings matrix (optional)

assay.used Name of assay used to generate DimReduc object

global Is this DimReduc global/persistent? If so, it will not be removed when removing its associated assay

stdev A vector of standard deviations

jackstraw A JackStrawData-class object associated with this DimReduc

misc A named list of unstructured miscellaneous data

key A one-length character vector with the object's key; keys must be one or more alphanumeric
```

characters followed by an underscore "_" (regex pattern "^[a-zA-Z][a-zA-Z0-9]*_\$")

See Also

Dimensional reduction object, validity, and interaction methods CreateDimReducObject(), DimReduc-validity, [.DimReduc(), [[.DimReduc(), dim.DimReduc(), merge.DimReduc(), print.DimReduc(), subset.DimReduc()

DimReduc-validity

Dimensional Reduction Validity

Description

Validation of DimReduc objects is handled by validObject

Cell Embeddings Validation

The cell embeddings matrix must be a numeric matrix of dimensions n_{cells} by $d_{dimensions}$; row names must be the cell names and column names must be the dimension identifier. The dimension identifier must be "key_dimension" (eg. "PC_1"). Dimension identifiers must be in order and cannot be skipped

Feature and Projected Feature Loadings Validation

blah

Standard Deviations Validation

blah

Distances 49

Key Validation

Keys must be a one-length character vector; a key must be composed of one of the following:

- An empty string (eg. "''") where nzchar() == 0
- An string composed of one or more alphanumeric values (both lower- and upper-case) that ends with an underscore ("_"); the first character must be a letter

Keys that are not empty strings are validated with the regex "^[a-zA-Z][a-zA-Z0-9]*_\$"

See Also

```
Dimensional reduction object, validity, and interaction methods CreateDimReducObject(), DimReduc-class, [.DimReduc(), [[.DimReduc(), dim.DimReduc(), merge.DimReduc(), print.DimReduc(), subset.DimReduc()
```

Distances

Get the Neighbor nearest neighbors distance matrix

Description

Get the Neighbor nearest neighbors distance matrix

Usage

```
Distances(object, ...)
## S3 method for class 'Neighbor'
Distances(object, ...)
```

Arguments

object An object

.. Arguments passed to other methods

Value

The distance matrix

50 droplevels.LogMap

droplevels.LogMap

Drop Unused Logical Map Values

Description

Remove any unused values from a logical map

Usage

```
## S3 method for class 'LogMap'
droplevels(x, ...)
```

Arguments

```
x A LogMap object
... Ignored
```

Value

x with values not present in any observation removed

See Also

```
Logical\ map\ objects,\ validity,\ and\ interaction\ methods:\ LogMap-validity,\ LogMap,\ as.\ matrix.\ LogMap(),\ intersect.\ LogMap(),\ labels.\ LogMap()
```

```
map <- LogMap(letters[1:10])
map[['obs']] <- c(1, 3, 7)
map[['entry']] <- c(2, 7, 10)

# Remove unused values
map <- droplevels(map)
map
map[[]]</pre>
```

Embeddings 51

Embeddings

Get Cell Embeddings

Description

Get Cell Embeddings

Usage

```
Embeddings(object, ...)
## S3 method for class 'DimReduc'
Embeddings(object, ...)
## S3 method for class 'Seurat'
Embeddings(object, reduction = "pca", ...)
```

Arguments

object An object

... Arguments passed to other methods

reduction Name of reduction to pull cell embeddings for

Value

The embeddings matrix

Examples

```
# Get the embeddings directly from a DimReduc object
Embeddings(object = pbmc_small[["pca"]])[1:5, 1:5]
# Get the embeddings from a specific DimReduc in a Seurat object
Embeddings(object = pbmc_small, reduction = "pca")[1:5, 1:5]
```

FetchData

Access cellular data

Description

Retrieves data (feature expression, PCA scores, metrics, etc.) for a set of cells in a Seurat object

52 FetchData

Usage

```
FetchData(object, ...)
## S3 method for class 'DimReduc'
FetchData(object, vars, cells = NULL, ...)
## S3 method for class 'Seurat'
FetchData(
   object,
   vars,
   cells = NULL,
   layer = NULL,
   clean = TRUE,
   slot = deprecated(),
   ...
)
```

Arguments

object	An object
	Arguments passed to other methods
vars	List of all variables to fetch, use keyword "ident" to pull identity classes
cells	Cells to collect data for (default is all cells)
layer	Layer to pull feature data for
clean	Remove cells that are missing data; choose from:
	• "all": consider all columns for cleaning
	 "ident": consider all columns except the identity class for cleaning
	• "project": consider all columns except the identity class for cleaning; fill missing identity values with the object's project
	• "none": do not clean
	Passing TRUE is a shortcut for "ident"; passing FALSE is a shortcut for "none"
slot	Deprecated in favor of layer

Value

A data frame with cells as rows and cellular data as columns

```
pc1 <- FetchData(object = pbmc_small, vars = 'PC_1')
head(x = pc1)
head(x = FetchData(object = pbmc_small, vars = c('groups', 'ident')))</pre>
```

FilterObjects 53

FilterObjects

Find Sub-objects of a Certain Class

Description

Get the names of objects within a Seurat object that are of a certain class

Usage

```
FilterObjects(object, classes.keep = c("Assay", "StdAssay", "DimReduc"))
```

Arguments

object A Seurat object

classes.keep A vector of names of classes to get

Value

A vector with the names of objects within the Seurat object that are of class classes. keep

Lifecycle

[Deprecated]

FilterObjects was deprecated in version 5.0.0; use .FilterObjects instead

Examples

FilterObjects(pbmc_small)

FOV-class

The Field of View Object

Description

A modern container for storing coordinates of spatially-resolved single cells. Capable of storing multiple cell segmentation boundary masks. Supports coordinates for spatially-resolved molecule (FISH) data. Compatible with SpatialImage

Slots

molecules A named list of Molecules objects defining spatially-resolved molecular coordinates boundaries A named list of Segmentation and Centroids objects defining spatially-resolved boundaries

assay A character naming the associated assay of the spatial coordinates

key A one-length character vector with the object's key; keys must be one or more alphanumeric characters followed by an underscore "_" (regex pattern "^[a-zA-Z][a-zA-Z0-9]*_\$")

54 FOV-methods

See Also

FOV-methods

FOV-methods

FOV Methods

Description

Methods for FOV objects

Usage

```
## S3 method for class 'FOV'
Cells(x, boundary = NULL, ...)
## S3 method for class 'FOV'
Features(x, set = NULL, ...)
## S3 method for class 'FOV'
FetchData(object, vars, cells = NULL, simplify = TRUE, ...)
## S3 method for class 'FOV'
GetTissueCoordinates(object, which = NULL, ...)
## S3 method for class 'FOV'
Keys(object, ...)
## S3 method for class 'FOV'
RenameCells(object, new.names = NULL, ...)
## S3 method for class 'FOV'
x$i, ...
## S3 method for class 'FOV'
x[i, j, ...]
## S3 method for class 'FOV'
x[[i, ...]]
## S3 method for class 'FOV'
length(x)
## S3 method for class 'FOV'
names(x)
## S3 method for class 'FOV'
subset(x, cells = NULL, features = NULL, ...)
```

FOV-methods 55

```
## S4 replacement method for signature 'FOV, character, missing, Centroids'
x[[i, j, ...]] <- value

## S4 replacement method for signature 'FOV, character, missing, Molecules'
x[[i, j, ...]] <- value

## S4 replacement method for signature 'FOV, character, missing, `NULL`'
x[[i, j, ...]] <- value

## S4 replacement method for signature 'FOV, character, missing, Segmentation'
x[[i, j, ...]] <- value

## S4 method for signature 'FOV'
show(object)</pre>
```

Arguments

x, object	A FOV object
boundary, set	Name of segmentation boundary or molecule set to extract cell or feature names for; pass NA to return all cells or feature names
• • •	Arguments passed to other methods
vars	A vector of variables to fetch; can be the name of a segmentation boundary, to get tissue coordinates, or molecule names, to get molecule coordinates
simplify	If only returning either boundary or molecule coordinates, return a single data frame instead of a list
which	Name of segmentation boundary or molecule set
new.names	vector of new cell names
i, cells	For [[and [[<-, the name of a segmentation or "molecules"; for FetchData, subset. and [, a vector of cells to keep
j, features	For subset and [, a vector of features to keep; for [[<-, not used
value	For [[<-, a replacement Molecules, Centroids, or Segmentation object; otherwise NULL to remove the boundary stored at i

Details

The following methods are defined for interacting with a FOV object:

Cells: Get cell names

Features: Get spatially-resolved molecule names

FetchData: Fetch boundary and/or molecule coordinates from a FOV object

GetTissueCoordinates: Get boundary or molecule coordinates from a FOV object

Keys: Get the keys of molecule sets contained within a FOV object

RenameCells: Update cell names

\$, [[: Extract a segmentation boundary

56 FOV-methods

length: Get the number of segmentation layers in a FOV object

names: Get the names of segmentation layers and molecule sets

subset, [: Subset a FOV object

[[<-: Add or remove segmentation layers and molecule information to/from a FOV object

show: Display an object summary to stdout

Value

Cells: A vector of cell names

Features: A vector of spatially-resolved molecule names; if no molecular information present, returns NULL

FetchData: If both molecule and boundary coordinates are requested, then a two-length list:

- "molecules": A data frame with the molecule coordinates requested. If molecules requested are keyed, the keys are preserved in the data frame
- "coordinates": A data frame with coordinates from the segmentation boundaries requested

If simplify is TRUE and only one data frame is generated, then only the data frame is returned. Otherwise, a one-length list is returned with the single data frame generated

GetTissueCoordinates: ...

Keys: A named vector of molecule set keys; names are the names of the molecule sets and values are the keys for the respective molecule set

RenameCells: object with the cells renamed to new.names

\$, [[: The segmentation boundary or spatially-resolved molecule information stored at i

length: The number of segmentation layers (Segmentation or Centroids objects)

names: A vector of segmentation boundary and molecule set names

subset: x with just the cells and features specified

[[<-: Varies depending on the class of value:

- If value is NULL, returns x with the boundary i removed; also allows removing molecules; does not allow removing the default segmentation
- If value is a Molecules, returns x with value stored in molecules; requires that i is "molecules"
- Otherwise, stores value as a segmentation boundary named i

show: Invisibly returns NULL

See Also

FOV-class

FOV-validity 57

FOV-validity

FOV Validity

Description

Validation of FOV objects is handled by validObject

Boundary Validation

blah

Molecule Validation

blah

See Also

validObject

GetImage

Get image data

Description

Get image data

Usage

```
GetImage(object, mode = c("grob", "raster", "plotly", "raw"), ...)
## S3 method for class 'Seurat'
GetImage(
  object,
  mode = c("grob", "raster", "plotly", "raw"),
  image = NULL,
  ...
)
```

Arguments

```
object An object

mode How to return the image; should accept one of "grob", "raster", "plotly", or "raw"

... Arguments passed to other methods

image Name of SpatialImage object to pull image data for; if NULL, will attempt to select an image automatically
```

58 GetTissueCoordinates

Value

```
Image data, varying depending on the value of mode:
```

```
"grob" An object representing image data inheriting from grob objects (eg. rastergrob)
```

"raster" An object of class raster

"plotly" A list with image data suitable for Plotly rendering, see plotly::layout for more details

"raw" The raw image data as stored in the object

See Also

layout

GetTissueCoordinates Get tissue coordinates

Description

Get tissue coordinates

Usage

```
GetTissueCoordinates(object, ...)
## S3 method for class 'Seurat'
GetTissueCoordinates(object, image = NULL, ...)
```

Arguments

object An object

... Arguments passed to other methods

image Name of SpatialImage object to get coordinates for; if NULL, will attempt to

select an image automatically

Value

A data frame with tissue coordinates

Graph-class 59

Graph-class

The Graph Class

Description

The Graph class inherits from dgCMatrix. We do this to enable future expandability of graphs.

Slots

assay.used Optional name of assay used to generate Graph object

See Also

```
dgCMatrix-class
Other graph: as.Graph()
```

HVFInfo

Highly Variable Features

Description

Get and set variable feature information for an Assay object. HVFInfo and VariableFeatures utilize generally variable features, while SVFInfo and SpatiallyVariableFeatures are restricted to spatially variable features

Usage

```
HVFInfo(object, method, status = FALSE, ...)
VariableFeatures(object, method = NULL, ...)
VariableFeatures(object, ...) <- value
SVFInfo(object, method, status, ...)
SpatiallyVariableFeatures(object, method, ...)
## S3 method for class 'Seurat'
HVFInfo(
   object,
   method = NULL,
   status = FALSE,
   assay = NULL,
   selection.method = deprecated(),</pre>
```

60 HVFInfo

```
)
## S3 method for class 'Seurat'
VariableFeatures(
  object,
 method = NULL,
  assay = NULL,
 nfeatures = NULL,
  layer = NA,
  simplify = TRUE,
  selection.method = deprecated(),
)
## S3 replacement method for class 'Seurat'
VariableFeatures(object, assay = NULL, ...) <- value</pre>
## S3 method for class 'Seurat'
SVFInfo(
  object,
 method = c("markvariogram", "moransi"),
  status = FALSE,
 assay = NULL,
  selection.method = deprecated(),
)
## S3 method for class 'Seurat'
SpatiallyVariableFeatures(
 object,
 method = "moransi",
  assay = NULL,
  decreasing = TRUE,
  selection.method = deprecated(),
)
## S3 method for class 'Assay'
HVFInfo(object, method, status = FALSE, selection.method = deprecated(), ...)
## S3 method for class 'Assay'
SpatiallyVariableFeatures(
 object,
 method = "moransi",
 decreasing = TRUE,
  selection.method = deprecated(),
)
```

HVFInfo 61

```
## S3 method for class 'Assay'
   SVFInfo(
      object,
     method = c("markvariogram", "moransi"),
     status = FALSE,
      selection.method = deprecated(),
   )
   ## S3 method for class 'Assay'
   VariableFeatures(object, method = NULL, selection.method = deprecated(), ...)
   ## S3 replacement method for class 'Assay'
   VariableFeatures(object, ...) <- value
   ## S3 method for class 'Assay5'
   HVFInfo(object, method = NULL, status = FALSE, layer = NULL, strip = TRUE, ...)
   ## S3 method for class 'Assay5'
   VariableFeatures(
     object,
     method = NULL,
     layer = NA,
      simplify = TRUE,
     nfeatures = Inf,
      selection.method = deprecated(),
   )
   ## S3 replacement method for class 'Assay5'
   VariableFeatures(object, method = "custom", layer = NULL, ...) <- value</pre>
Arguments
   object
                    An object
   method
                    Which method to pull. For HVFInfo and VariableFeatures, choose one from
                    one of the following:
                      • "vst"
                      • "sctransform" or "sct"
                      • "mean.var.plot", "dispersion", "mvp", or "disp"
                    For SVFInfo and Spatially Variable Features, choose from:
                      · "markvariogram"
                      · "moransi"
                    Add variable status to the resulting data frame
   status
                    Arguments passed to other methods
```

62 Idents

value A character vector of variable features

assay Name of assay to pull highly variable feature information for

selection.method

[Deprecated]

nfeatures Maximum number of features to select when simplifying

layer Layer to pull variable features for

simplify When pulling for multiple layers, combine into a single vector and select a com-

mon set of variable features for all layers

decreasing Return features in decreasing order (most spatially variable first).

strip Remove method/layer identifiers from highly variable data frame

Value

HVFInfo: A data frame with feature means, dispersion, and scaled dispersion

VariableFeatures: a vector of the variable features

SVFInfo: a data frame with the spatially variable features

SpatiallyVariableFeatures: a character vector of the spatially variable features

Examples

```
# Get the HVF info from a specific Assay in a Seurat object
HVFInfo(object = pbmc_small, assay = "RNA")[1:5, ]
# Get the HVF info directly from an Assay object
HVFInfo(pbmc_small[["RNA"]], method = 'vst')[1:5, ]
```

Idents

Get, set, and manipulate an object's identity classes

Description

Get, set, and manipulate an object's identity classes

Usage

```
Idents(object, ...)
Idents(object, ...) <- value
RenameIdents(object, ...)
ReorderIdent(object, var, ...)
SetIdent(object, ...)</pre>
```

Idents 63

```
StashIdent(object, save.name, ...)
## S3 method for class 'Seurat'
Idents(object, ...)
## S3 replacement method for class 'Seurat'
Idents(object, cells = NULL, drop = FALSE, replace = FALSE, ...) <- value</pre>
## S3 method for class 'Seurat'
ReorderIdent(
  object,
 var,
  reverse = FALSE,
  afxn = mean,
  reorder.numeric = FALSE,
)
## S3 method for class 'Seurat'
RenameIdents(object, ...)
## S3 method for class 'Seurat'
SetIdent(object, cells = NULL, value, ...)
## S3 method for class 'Seurat'
StashIdent(object, save.name = "orig.ident", ...)
## S3 method for class 'Seurat'
droplevels(x, ...)
## S3 method for class 'Seurat'
levels(x)
## S3 replacement method for class 'Seurat'
levels(x) <- value
```

Arguments

• • •	Arguments passed to other methods; for RenameIdents: named arguments as old.ident = new.ident; for ReorderIdent: arguments passed on to FetchData
value	The name of the identities to pull from object metadata or the identities themselves
var	Feature or variable to order on
save.name	Store current identity information under this name
cells	Set cell identities for specific cells
drop	Drop unused levels

64 Idents

replace Replace identities for unset cells with NA

reverse Reverse ordering

afxn Function to evaluate each identity class based on; default is mean

reorder.numeric

Rename all identity classes to be increasing numbers starting from 1 (default is

FALSE)

x, object An object

Value

Idents: The cell identities

Idents<-: object with the cell identities changed

RenameIdents: An object with selected identity classes renamed

ReorderIdent: An object with

SetIdent: An object with new identity classes set StashIdent: An object with the identities stashed

```
# Get cell identity classes
Idents(pbmc_small)
# Set cell identity classes
# Can be used to set identities for specific cells to a new level
Idents(pbmc_small, cells = 1:4) <- 'a'</pre>
head(Idents(pbmc_small))
# Can also set idents from a value in object metadata
colnames(pbmc_small[[]])
Idents(pbmc_small) <- 'RNA_snn_res.1'</pre>
levels(pbmc_small)
# Rename cell identity classes
# Can provide an arbitrary amount of idents to rename
levels(pbmc_small)
pbmc_small <- RenameIdents(pbmc_small, '0' = 'A', '2' = 'C')</pre>
levels(pbmc_small)
## Not run:
head(Idents(pbmc_small))
pbmc_small <- ReorderIdent(pbmc_small, var = 'PC_1')</pre>
head(Idents(pbmc_small))
## End(Not run)
# Set cell identity classes using SetIdent
cells.use <- WhichCells(pbmc_small, idents = '1')</pre>
pbmc_small <- SetIdent(pbmc_small, cells = cells.use, value = 'B')</pre>
```

Images 65

```
head(pbmc_small[[]])
pbmc_small <- StashIdent(pbmc_small, save.name = 'idents')
head(pbmc_small[[]])

# Get the levels of identity classes of a Seurat object
levels(x = pbmc_small)

# Reorder identity classes
levels(x = pbmc_small)
levels(x = pbmc_small) <- c('C', 'A', 'B')
levels(x = pbmc_small)</pre>
```

Images

Pull spatial image names

Description

List the names of SpatialImage objects present in a Seurat object. If assay is provided, limits search to images associated with that assay

Usage

```
Images(object, assay = NULL)
```

Arguments

object A Seurat object

assay Name of assay to limit search to

Value

A list of image names

```
## Not run:
Images(object)
## End(Not run)
```

66 Indices

Index

Get Neighbor algorithm index

Description

Get Neighbor algorithm index

Usage

```
Index(object, ...)
Index(object, ...) <- value
## S3 method for class 'Neighbor'
Index(object, ...)
## S3 replacement method for class 'Neighbor'
Index(object, ...) <- value</pre>
```

Arguments

object An object

... Arguments passed to other methods

value The index to store

Value

Returns the value in the alg.idx slot of the Neighbor object

Idents<-: A Neighbor object with the index stored

Indices

Get Neighbor nearest neighbor index matrices

Description

Get Neighbor nearest neighbor index matrices

Usage

```
Indices(object, ...)
## S3 method for class 'Neighbor'
Indices(object, ...)
```

intersect.LogMap 67

Arguments

object An object
... Arguments passed to other methods

Value

A matrix with the nearest neighbor indices

intersect.LogMap

Find Common Logical Map Values

Description

Identify values in a logical map that are common to every observation

Usage

```
## S3 method for class 'LogMap'
intersect(x, y = missing_arg(), ...)
```

Arguments

```
x A LogMap objecty Ignored... Ignored
```

Value

The values of x that are present in **every** observation

See Also

```
Logical map objects, validity, and interaction methods: LogMap-validity, LogMap, as.matrix.LogMap(), droplevels.LogMap(), labels.LogMap()
```

```
map <- LogMap(letters[1:10])
map[['obs']] <- c(1, 3, 7)
map[['entry']] <- c(2, 7, 10)

# Identify values that are present in every observation
intersect(map)</pre>
```

68 IsMatrixEmpty

IsGlobal

Is an object global/persistent?

Description

Typically, when removing Assay objects from an Seurat object, all associated objects (eg. DimReduc, Graph, and SeuratCommand objects) are removed as well. If an associated object is marked as global/persistent, the associated object will remain even if its original assay was deleted

Usage

```
IsGlobal(object, ...)
## Default S3 method:
IsGlobal(object, ...)
## S3 method for class 'DimReduc'
IsGlobal(object, ...)
```

Arguments

object An object

... Arguments passed to other methods

Value

TRUE if the object is global/persistent otherwise FALSE

Examples

```
IsGlobal(pbmc_small[['pca']])
```

IsMatrixEmpty

Check if a matrix is empty

Description

Takes a matrix and asks if it's empty (either 0x0 or 1x1 with a value of NA)

Usage

```
IsMatrixEmpty(x)
## Default S3 method:
IsMatrixEmpty(x)
```

IsNamedList 69

Arguments

x A matrix

Value

Whether or not x is empty

Examples

```
IsMatrixEmpty(new("matrix"))
IsMatrixEmpty(matrix())
IsMatrixEmpty(matrix(1:3))
```

IsNamedList

Check List Names

Description

Check to see if a list has names; also check to enforce that all names are present and unique

Usage

```
IsNamedList(x, all.unique = TRUE, allow.empty = FALSE, pass.zero = FALSE)
```

Arguments

```
x A list

all.unique Require that all names are unique from one another

allow.empty Allow empty (nchar = 0) names

pass.zero Pass on zero-length lists
```

Value

```
TRUE if ..., otherwise FALSE
```

```
IsNamedList(list())
IsNamedList(list(), pass.zero = TRUE)
IsNamedList(list(1, 2, 3))
IsNamedList(list(a = 1, b = 2, c = 3))
IsNamedList(list(a = 1, 2, c = 3))
IsNamedList(list(a = 1, 2, c = 3), allow.empty = TRUE)
IsNamedList(list(a = 1, a = 2, a = 3))
IsNamedList(list(a = 1, a = 2, a = 3), all.unique = FALSE)
```

70 JackStrawData-methods

JackStrawData-class The JackStrawData Class

Description

The JackStrawData is used to store the results of a JackStraw computation.

Slots

```
empirical.p.values Empirical p-values
fake.reduction.scores Fake reduction scores
empirical.p.values.full Empirical p-values on full
overall.p.values Overall p-values from ScoreJackStraw
```

JackStrawData-methods JackStrawData Methods

Description

Methods for JackStrawData objects for generics defined in other packages

Usage

```
## S3 method for class 'JackStrawData'
.DollarNames(x, pattern = "")
## S3 method for class 'JackStrawData'
x$i, ...
## S3 method for class 'JackStrawData'
as.logical(x, ...)
## S4 method for signature 'JackStrawData'
show(object)
```

Arguments

```
    x, object A JackStrawData object
    pattern A regular expression. Only matching names are returned.
    i A JackStrawData slot name
    ... Ignored
```

JoinLayers 71

Value

```
$: Slot i from x as.logical: TRUE if empirical p-values have been calculated otherwise FALSE show: Prints summary to stdout and invisibly returns NULL
```

Functions

- .DollarNames(JackStrawData): Autocompletion for \$ access on a JackStrawData object
- \$: Access data from a JackStrawData object
- as.logical(JackStrawData): Have empirical p-values for a JackStrawData object been calculated
- show(JackStrawData): Overview of a JackStrawData object

JoinLayers

Split and Join Layers Together

Description

Split and Join Layers Together

Usage

```
JoinLayers(object, ...)
## S3 method for class 'Assay5'
JoinLayers(object, layers = NULL, new = NULL, ...)
## S3 method for class 'Seurat'
JoinLayers(object, assay = NULL, layers = NULL, new = NULL, ...)
```

Arguments

object An object
... Arguments passed to other methods
layers Names of layers to split or join
new Name of new layers
assay Name of assay to split layers

Value

object with the layers specified joined

Get and set JackStraw information

JS

Description

Get and set JackStraw information

Usage

```
JS(object, ...)
JS(object, ...) <- value

## S3 method for class 'JackStrawData'
JS(object, slot, ...)

## S3 replacement method for class 'JackStrawData'
JS(object, slot, ...) <- value

## S3 method for class 'DimReduc'
JS(object, slot = NULL, ...)

## S3 replacement method for class 'DimReduc'
JS(object, slot = NULL, ...) <- value</pre>
```

Arguments

object	An object
	Arguments passed to other methods
value	JackStraw information
slot	Name of slot to store JackStraw scores to Can shorten to 'empirical', 'fake', 'full', or 'overall'

Value

```
JS: either a JackStrawData object or the specified jackstraw data JS<-: object with the update jackstraw information
```

Key 73

Key

Get and set object keys

Description

Get and set object keys

Usage

```
Key(object, ...)
Keys(object, ...)
Key(object, ...) <- value</pre>
## S3 method for class 'Assay'
Key(object, ...)
## S3 replacement method for class 'Assay'
Key(object, ...) <- value</pre>
## S3 method for class 'Assay5'
Key(object, ...)
## S3 replacement method for class 'Assay5'
Key(object, ...) <- value</pre>
## S3 method for class 'DimReduc'
Key(object, ...)
## S3 replacement method for class 'DimReduc'
Key(object, ...) <- value</pre>
## S3 method for class 'Seurat'
Key(object, ...)
## S3 method for class 'Seurat'
Keys(object, ...)
```

Arguments

object An object
... Arguments passed to other methods
value Key value

74 labels.LogMap

Value

```
Key: the object key

Keys: a named vector of keys of sub-objects

Key<-: object with an updated key
```

Examples

```
# Get an Assay key
Key(pbmc_small[["RNA"]])

# Set the key for an Assay
Key(pbmc_small[["RNA"]]) <- "newkey_"
Key(pbmc_small[["RNA"]])

# Get a DimReduc key
Key(object = pbmc_small[["pca"]])

# Set the key for DimReduc
Key(object = pbmc_small[["pca"]]) <- "newkey2_"
Key(object = pbmc_small[["pca"]])

# Show all keys associated with a Seurat object
Key(object = pbmc_small)
Keys(object = pbmc_small)</pre>
```

labels.LogMap

Find Observations by Value

Description

Identify the observations that contain a specific value in a logical map

```
## S3 method for class 'LogMap'
labels(
  object,
  values,
  select = c("first", "last", "common", "all"),
  simplify = TRUE,
  ...
)
```

LayerData 75

Arguments

Value

labels: A list, or vector if simplify is TRUE, of all values and the observations they're found in, according to the value of select

See Also

Logical map objects, validity, and interaction methods: LogMap-validity, LogMap, as.matrix.LogMap(), droplevels.LogMap(), intersect.LogMap()

Examples

```
map <- LogMap(letters[1:10])
map[['obs']] <- c(1, 3, 7)
map[['entry']] <- c(2, 7, 10)

# Find observations for a set of values labels(map, c('a', 'b', 'g'))</pre>
```

LayerData

Query and Manipulate Assay Layers

Description

Query and Manipulate Assay Layers

76 LayerData

```
LayerData(object, layer, ...)
LayerData(object, layer, ...) <- value
Layers(object, ...)
## S3 method for class 'Assay'
LayerData(
 object,
  layer = NULL,
  cells = NULL,
  features = NULL,
  slot = deprecated(),
)
## S3 replacement method for class 'Assay'
LayerData(object, layer, ...) <- value
## S3 method for class 'Assay'
Layers(object, search = NA, ...)
## S3 method for class 'Assay5'
LayerData(
  object,
  layer = NULL,
  cells = NULL,
  features = NULL,
  fast = FALSE,
  slot = deprecated(),
  . . .
)
## S3 replacement method for class 'Assay5'
LayerData(object, layer, features = NULL, cells = NULL, ...) <- value
## S3 method for class 'Assay5'
Layers(object, search = NA, ...)
## S3 method for class 'Seurat'
LayerData(object, layer = NULL, assay = NULL, slot = deprecated(), ...)
## S3 replacement method for class 'Seurat'
LayerData(object, layer, assay = NULL, ...) <- value
## S3 method for class 'Seurat'
Layers(object, search = NA, assay = NULL, ...)
```

Loadings 77

Arguments

object An object

layer Name of layer to fetch or set

... Arguments passed to other methods

value New two-dimensional data to be added as a layer

features, cells

Vectors of features/cells to include

slot [Deprecated]

search A pattern to search layer names for; pass one of:

• "NA" to pull all layers

• "NULL" to pull the default layer(s)

• a regular expression that matches layer names

fast Determine how to return the layer data; choose from:

FALSE Apply any transpositions and attempt to add feature/cell names (if sup-

ported) back to the layer data

NA Attempt to add feature/cell names back to the layer data, skip any transposi-

tions

TRUE Do not apply any transpositions or add feature/cell names to the layer data

assay Name of assay to fetch layer data from or assign layer data to

Value

LayerData: the layer data for layer from object

Layer<-: object with value added as a layer named layer

Layers: the names of the layers present in object

Loadings

Get and set feature loadings

Description

Get and set feature loadings

```
Loadings(object, ...)
Loadings(object, ...) <- value
## S3 method for class 'DimReduc'
Loadings(object, projected = FALSE, ...)</pre>
```

78 LogMap

```
## S3 replacement method for class 'DimReduc'
Loadings(object, projected = TRUE, ...) <- value
## S3 method for class 'Seurat'
Loadings(object, reduction = "pca", projected = FALSE, ...)</pre>
```

Arguments

object An object

... Arguments passed to other methods

value Feature loadings to add

projected Pull the projected feature loadings?

reduction Name of reduction to pull feature loadings for

Value

```
Loadings: the feature loadings for object
Loadings<-: object with the updated loadings
```

Examples

```
# Get the feature loadings for a given DimReduc
Loadings(object = pbmc_small[["pca"]])[1:5,1:5]

# Set the feature loadings for a given DimReduc
new.loadings <- Loadings(object = pbmc_small[["pca"]])
new.loadings <- new.loadings + 0.01
Loadings(object = pbmc_small[["pca"]]) <- new.loadings
# Get the feature loadings for a specified DimReduc in a Seurat object
Loadings(object = pbmc_small, reduction = "pca")[1:5,1:5]</pre>
```

LogMap

A Logical Map

Description

A simple container for storing mappings of values using logical matrices. Keeps track of which values (rows) are present in which observations (columns). LogMap objects can be created with LogMap(); queries can be performed with [[and observations can be added or removed with [[<-

LogMap 79

Usage

```
LogMap(y)

## S4 method for signature 'LogMap,character,missing'
x[[i, j, ...]]

## S4 method for signature 'LogMap,missing,missing'
x[[i, j, ...]]

## S4 method for signature 'LogMap,`NULL`,missing'
x[[i, j, ...]]

## S4 replacement method for signature 'LogMap,character,missing,character'
x[[i, j, ...]] <- value

## S4 replacement method for signature 'LogMap,character,missing,integer'
x[[i, j, ...]] <- value

## S4 replacement method for signature 'LogMap,character,missing,`NULL`'
x[[i, j, ...]] <- value

## S4 replacement method for signature 'LogMap,character,missing,`NULL`'
x[[i, j, ...]] <- value
```

Arguments

У	A character vector
X	A LogMap object
i	A character vector of length 1, or NULL
j	Not used
	Ignored
value	A character or integer vector of values to record in the map for i , or NULL to remove the record for i

Value

LogMap: A new LogMap object with zero columns and length(x = x) rows; rownames are set to x [[: if i is a character vector, the rownames that are mapped to i; otherwise the rownames of x [[<-: If value is NULL, then x without the observations for i; otherwise, x with a new column for i recording a TRUE for all values present in value

Slots

.Data A logical matrix with at least one row

80 LogMap-validity

See Also

Logical map objects, validity, and interaction methods: LogMap-validity, as.matrix.LogMap(), droplevels.LogMap(), intersect.LogMap(), labels.LogMap()

Examples

```
# Create a LogMap
map <- LogMap(letters[1:10])</pre>
# Get the names of values in the LogMap
map[[NULL]]
rownames(map)
# Add an observation to the LogMap
map[['obs']] \leftarrow c(1, 3, 7)
map[['entry']] \leftarrow c(2, 7, 10)
map
# Get the names of observations in the LogMap
colnames(map)
# Fetch an observation from the LogMap
map[['obs']]
# Get the full logical matrix
map[[]]
# Remove an observation from the LogMap
map[['obs']] <- NULL</pre>
map[['entry']] <- NULL</pre>
map
```

LogMap-validity

Logical Map Validity

Description

Validation of LogMap objects is handled by validObject

Data Validation

Logical maps must be a logical matrix containing only TRUE or FALSE values

Value Validation

All values must be named within the rownames of the object. Duplicate or empty ("") values are not allowed

LogSeuratCommand 81

Observation Validation

All observations must be named within the column names of the object. Duplicate or empty ("") observations are not allowed

See Also

```
validObject
```

Logical map objects, validity, and interaction methods: LogMap, as.matrix.LogMap(), droplevels.LogMap(), intersect.LogMap(), labels.LogMap()

Examples

```
map <- LogMap(letters[1:10])
map[['obs']] <- c(1, 3, 7)
map[['entry']] <- c(2, 7, 10)
validObject(map)</pre>
```

LogSeuratCommand

Log a command

Description

Logs command run, storing the name, timestamp, and argument list. Stores in the Seurat object

Usage

```
LogSeuratCommand(object, return.command = FALSE)
```

Arguments

```
object Name of Seurat object
return.command Return a SeuratCommand object instead
```

Value

If return.command, returns a SeuratCommand object; otherwise, returns the Seurat object with command stored

See Also

Command

```
Command log object and interaction methods $.SeuratCommand(), .DollarNames.SeuratCommand(), SeuratCommand-class, [.SeuratCommand(), as.list.SeuratCommand()
```

82 merge.Assay

merge.Assay

Merge Assays

Description

Merge one or more v3 assays together

Usage

```
## S3 method for class 'Assay'
merge(
  x = NULL,
  y = NULL,
  add.cell.ids = NULL,
  merge.data = TRUE,
  labels = NULL,
  collapse = TRUE,
  ...
)
```

Arguments

```
x An Assay object
y One or more Assay objects
add.cell.ids A character vector of length(x = c(x, y)); appends the corresponding values to the start of each objects' cell names

merge.data Merge the data slots instead of just merging the counts (which requires renormalization); this is recommended if the same normalization approach was applied to all objects

labels, collapse

Currently unused
... Ignored
```

Value

A new assay with data merged from c(x, y)

See Also

```
v3 Assay object, validity, and interaction methods: $.Assay(), Assay-class, Assay-validity, CreateAssayObject(), [.Assay(), [.Assay(), dim.Assay(), dimnames.Assay(), split.Assay(), subset.Assay()
```

merge.Assay5

Description

Merge one or more v5 assays together

Usage

```
## S3 method for class 'Assay5'
merge(x, y, labels = NULL, add.cell.ids = NULL, collapse = FALSE, ...)
```

Arguments

X	An Assay5 object
У	One or more Assay5 objects
labels	A character vector equal to the number of objects; defaults to as.character(seq_along(c(x, y)))
add.cell.ids	A character vector equal to the number of objects provided to append to all cell names; if TRUE, uses labels as add.cell.ids
collapse	If TRUE, merge layers of the same name together; if FALSE, appends labels to the layer name
	Ignored

Details

Note: collapsing layers is currently not supported

Value

A new v5 assay with data merged from c(x, y)

See Also

```
v5 Assay object, validity, and interaction methods: $.Assay5(), Assay5-class, Assay5-validity, [.Assay5(), [[.Assay5(), dim.Assay5(), dimnames.Assay5(), split.Assay5(), subset.Assay5()
```

84 merge.Seurat

merge.DimReduc

Merge Dimensional Reductions

Description

Merge two or more dimensional reduction together

Usage

```
## S3 method for class 'DimReduc'
merge(x = NULL, y = NULL, add.cell.ids = NULL, ...)
```

Arguments

```
    x A DimReduc object
    y One or more DimReduc objects
    add.cell.ids A character vector equal to the number of objects provided to append to all cell names; if TRUE, uses labels as add.cell.ids
    ... Ignored
```

Value

A new DimReduc object with data merged from c(x, y)

See Also

Dimensional reduction object, validity, and interaction methods CreateDimReducObject(), DimReduc-class, DimReduc-validity, [.DimReduc(), [[.DimReduc(), dim.DimReduc(), print.DimReduc(), subset.DimReduc()

merge.Seurat

Merge Seurat Objects

Description

Merge Seurat Objects

```
## S3 method for class 'Seurat'
merge(
  x = NULL,
  y = NULL,
  add.cell.ids = NULL,
  collapse = FALSE,
  merge.data = TRUE,
```

merge.Seurat 85

```
merge.dr = FALSE,
project = getOption(x = "Seurat.object.project", default = "SeuratProject"),
...
)
```

Arguments

X	A Seurat object
У	A single Seurat object or a list of Seurat objects
add.cell.ids	A character vector of length($x = c(x, y)$); appends the corresponding values to the start of each objects' cell names
collapse	If TRUE, merge layers of the same name together; if FALSE, appends labels to the layer name $$
merge.data	Merge the data slots instead of just merging the counts (which requires renormalization); this is recommended if the same normalization approach was applied to all objects
merge.dr	Choose how to handle merging dimensional reductions:
	• "TRUE": merge dimensional reductions with the same name across objects; dimensional reductions with different names are added as-is
	 "NA": keep dimensional reductions from separate objects separate; will append the project name for duplicate reduction names
	"FALSE": do not add dimensional reductions
project	Project name for the Seurat object
	Arguments passed to other methods

Value

merge: Merged object

Merge Details

When merging Seurat objects, the merge procedure will merge the Assay level counts and potentially the data slots (depending on the merge.data parameter). It will also merge the cell-level meta data that was stored with each object and preserve the cell identities that were active in the objects pre-merge. The merge will optionally merge reductions depending on the values passed to merge.dr if they have the same name across objects. Here the embeddings slots will be merged and if there are differing numbers of dimensions across objects, only the first N shared dimensions will be merged. The feature loadings slots will be filled by the values present in the first object. The merge will not preserve graphs, logged commands, or feature-level metadata that were present in the original objects. If add.cell.ids isn't specified and any cell names are duplicated, cell names will be appended with X, where X is the numeric index of the object in C(X, Y).

See Also

```
Seurat object, validity, and interaction methods $.Seurat(), Seurat-class, Seurat-validity,
[[.Seurat(), [[<-,Seurat, NULL, [[<-,Seurat, dim.Seurat(), dimnames.Seurat(), names.Seurat(),
subset.Seurat()</pre>
```

86 Misc

Examples

```
# `merge' examples
# merge two objects
merge(pbmc_small, y = pbmc_small)
# to merge more than two objects, pass one to x and a list of objects to y
merge(pbmc_small, y = c(pbmc_small, pbmc_small))
```

Misc

Get and set miscellaneous data

Description

Get and set miscellaneous data

```
Misc(object, ...)
Misc(object, ...) <- value
## S3 method for class 'Assay'
Misc(object, slot = NULL, ...)
## S3 replacement method for class 'Assay'
Misc(object, slot, ...) <- value
## S3 method for class 'Assay5'
Misc(object, slot = NULL, ...)
## S3 replacement method for class 'Assay5'
Misc(object, slot, ...) <- value
## S3 method for class 'DimReduc'
Misc(object, slot = NULL, ...)
## S3 replacement method for class 'DimReduc'
Misc(object, slot, ...) <- value
## S3 method for class 'Seurat'
Misc(object, slot = NULL, ...)
## S3 replacement method for class 'Seurat'
Misc(object, slot, ...) <- value
```

Molecules-class 87

Arguments

object An object

... Arguments passed to other methods

value Data to add

slot Name of specific bit of meta data to pull

Value

Miscellaneous data

An object with miscellaneous data added

Examples

```
# Get the misc info
Misc(object = pbmc_small, slot = "example")
# Add misc info
Misc(object = pbmc_small, slot = "example") <- "testing_misc"</pre>
```

Molecules-class

The Spatial Molecules Class

Description

The Spatial Molecules Class

Slots

```
.Data A list of SpatialPoints objects key The key for the Molecules
```

See Also

Molecules methods: Molecules-methods

 $Segmentation\ layer\ classes:\ Centroids-class,\ Centroids-methods,\ Molecules-methods,\ Segmentation-class,\ Segmentation-methods$

88 Molecules-methods

Molecules-methods

Molecules Methods

Description

Methods for Molecules objects

Usage

```
## S3 method for class 'Molecules'
Features(x, ...)

## S3 method for class 'Molecules'
GetTissueCoordinates(object, features = NULL, ...)

## S3 method for class 'Molecules'
subset(x, features = NULL, ...)

## S4 method for signature 'Molecules'
show(object)
```

Arguments

x, object A Molecules object

... Arguments passed to other methods

features A vector of molecule names to keep; if NULL, defaults to all molecules

Details

Features: Get spatially-resolved molecule names

GetTissueCoordinates: Get spatially-resolved molecule coordinates

subset: Subset a Molecules object to certain molecules

show: Display an object summary to stdout

Value

Features: A vector of spatially-resolved molecule names; if no molecular information present, returns NULL

GetTissueCoordinates: A data frame with three columns:

- "x": the x-coordinate of a molecule
- "y": the y-coordinate of a molecule
- "molecule": the molecule name

subset: x subsetted to the features specified by features

show: Invisibly returns NULL

names.Seurat 89

See Also

```
Molecules-class
```

Segmentation layer classes: Centroids-class, Centroids-methods, Molecules-class, Segmentation-class, Segmentation-methods

names.Seurat

Subobject Names

Description

Get the names of subobjects within a Seurat object

Usage

```
## S3 method for class 'Seurat'
names(x)
```

Arguments

Χ

A Seurat object

Value

The names of all of the following subobjects within x:

- v3 and v5 assays
- dimensional reductions
- images and FOVs
- nearest-neighbor graphs

See Also

```
Seurat object, validity, and interaction methods $.Seurat(), Seurat-class, Seurat-validity,
[[.Seurat(), [[<-,Seurat, NULL, [[<-,Seurat, dim.Seurat(), dimnames.Seurat(), merge.Seurat(),
subset.Seurat()</pre>
```

Examples

```
names(pbmc_small)
```

90 Neighbor-methods

Neighbor-class

The Neighbor class

Description

The Neighbor class is used to store the results of neighbor finding algorithms

Slots

```
nn.idx Matrix containing the nearest neighbor indices
```

nn.dist Matrix containing the nearest neighbor distances

alg.idx The neighbor finding index (if applicable). E.g. the annoy index

alg.info Any information associated with the algorithm that may be needed downstream (e.g. distance metric used with annoy is needed when reading in from stored file).

cell. names Names of the cells for which the neighbors have been computed.

Neighbor-methods

Neighbor Methods

Description

Methods for Neighbor objects for generics defined in other packages

Usage

```
## S3 method for class 'Neighbor'
dim(x)
## S4 method for signature 'Neighbor'
show(object)
```

Arguments

```
x, object A Neighbor object
```

Value

dim Dimensions of the indices matrix

show: Prints summary to stdout and invisibly returns NULL

Functions

- dim(Neighbor): Dimensions of the neighbor indices
- show(Neighbor): Overview of a Neighbor object

Overlay 91

Overlay

Overlay Spatial Objects Over One Another

Description

Create an overlay of some query spatial object (x) against some target object (y). Basically, find all components of a query that fall within the bounds of a target spatial region

Usage

```
Overlay(x, y, invert = FALSE, ...)
## S4 method for signature 'Centroids,SpatialPolygons'
Overlay(x, y, invert = FALSE, ...)
## S4 method for signature 'Segmentation,SpatialPolygons'
Overlay(x, y, invert = FALSE, ...)
## S4 method for signature 'Molecules,SpatialPolygons'
Overlay(x, y, invert = FALSE, ...)
## S4 method for signature 'FOV,Spatial'
Overlay(x, y, invert = FALSE, ...)
## S4 method for signature 'FOV,SpatialPolygons'
Overlay(x, y, invert = FALSE, ...)
## S4 method for signature 'FOV,FOV'
Overlay(x, y, invert = FALSE, ...)
```

Arguments

X	Query Spatial object
У	Target Spatial object
invert	Invert the overlay and return only the components of x that fall $outside$ the bounds of y
	Ignored

Value

x with only the components that fall within the bounds of y

Note

This function requires the sf package to be installed

92 pbmc_small

PackageCheck

Check the existence of a package

Description

Check the existence of a package

Usage

```
PackageCheck(..., error = TRUE)
```

Arguments

.. Package names

error If true, throw an error if the package doesn't exist

Value

Invisibly returns boolean denoting if the package is installed

Lifecycle

[Deprecated]

PackageCheck was deprecated in version 5.0.0; please use rlang::check_installed() instead

Examples

```
PackageCheck("SeuratObject", error = FALSE)
```

pbmc_small

A small example version of the PBMC dataset

Description

A subsetted version of 10X Genomics' 3k PBMC dataset

```
pbmc_small
```

print.DimReduc 93

Format

A Seurat object with the following slots filled

assays Currently only contains one assay ("RNA" - scRNA-seq expression data)

counts - Raw expression data

- data Normalized expression data
- scale.data Scaled expression data
- var.features names of the current features selected as variable
- meta.features Assay level metadata such as mean and variance

meta.data Cell level metadata

active.assay Current default assay

active.ident Current default idents

graphs Neighbor graphs computed, currently stores the SNN

reductions Dimensional reductions: currently PCA and tSNE

version Seurat version used to create the object

commands Command history

Source

https://support.10xgenomics.com/single-cell-gene-expression/datasets/1.1.0/pbmc3k

print.DimReduc

Print Top Feature Loadings

Description

Prints a set of features that most strongly define a set of components; **note**: requires feature loadings to be present in order to work

Usage

```
## S3 method for class 'DimReduc'
print(x, dims = 1:5, nfeatures = 20, projected = FALSE, ...)
```

Arguments

x A DimReduc object

dims Number of dimensions to display

nfeatures Number of genes to display

projected Use projected slot

... Ignored

94 Project

Value

Displays set of features defining the components and invisibly returns x

See Also

cat

Dimensional reduction object, validity, and interaction methods CreateDimReducObject(), DimReduc-class, DimReduc-validity, [.DimReduc(), [[.DimReduc(), dim.DimReduc(), merge.DimReduc(), subset.DimReduc()

Examples

```
pca <- pbmc_small[["pca"]]
print(pca)</pre>
```

Project

Get and set project information

Description

Get and set project information

Usage

```
Project(object, ...)
Project(object, ...) <- value

## S3 method for class 'Seurat'
Project(object, ...)

## S3 replacement method for class 'Seurat'
Project(object, ...) <- value</pre>
```

Arguments

object An object

... Arguments passed to other methods

value Project information to set

Value

Project information

An object with project information added

Radius 95

Radius

Get the spot radius from an image

Description

Get the spot radius from an image

Usage

```
Radius(object)
```

Arguments

object

An image object

Value

The radius size

RandomName

Generate a random name

Description

Make a name from randomly sampled characters, pasted together with no spaces

Usage

```
RandomName(length = 5L, chars = letters, ...)
```

Arguments

length How long should the name be

chars A vector of 1-length characters to use to generate the name

... Extra parameters passed to sample

Value

A character with nchar == length of randomly sampled letters

See Also

sample

96 RenameAssays

Examples

```
set.seed(42L)
RandomName()
RandomName(7L, replace = TRUE)
```

RenameAssays

Rename assays in a Seurat object

Description

Rename assays in a Seurat object

Usage

```
RenameAssays(
  object,
  assay.name = NULL,
  new.assay.name = NULL,
  verbose = TRUE,
  ...
)
```

Arguments

```
object A Seurat object
assay.name original name of assay
new.assay.name new name of assay
verbose Whether to print messages
... Named arguments as old.assay = new.assay
```

Value

object with assays renamed

Examples

```
RenameAssays(object = pbmc_small, RNA = 'rna')
```

RenameCells 97

RenameCells Rename cells

Description

Change the cell names in all the different parts of an object. Can be useful before combining multiple objects.

Usage

```
RenameCells(object, ...)
## S3 method for class 'Assay'
RenameCells(object, new.names = NULL, ...)
## S3 method for class 'Assay5'
RenameCells(object, new.names = NULL, ...)
## S3 method for class 'DimReduc'
RenameCells(object, new.names = NULL, ...)
## S3 method for class 'Neighbor'
RenameCells(object, old.names = NULL, new.names = NULL, ...)
## S3 method for class 'Seurat'
RenameCells(
 object,
  add.cell.id = missing_arg(),
 new.names = missing_arg(),
  for.merge = deprecated(),
)
```

Arguments

object	An object
	Arguments passed to other methods
new.names	vector of new cell names
old.names	vector of old cell names
add.cell.id	prefix to add cell names
for.merge	Deprecated

Details

If add.cell.id is set a prefix is added to existing cell names. If new.names is set these will be used to replace existing names.

Value

An object with new cell names

Examples

```
# Rename cells in an Assay
head(x = colnames(x = pbmc_small[["RNA"]]))
renamed.assay <- RenameCells(</pre>
    pbmc_small[["RNA"]],
    new.names = paste0("A_", colnames(x = pbmc_small[["RNA"]]))
head(x = colnames(x = renamed.assay))
# Rename cells in a DimReduc
head(x = Cells(x = pbmc_small[["pca"]]))
renamed.dimreduc <- RenameCells(</pre>
    object = pbmc_small[["pca"]],
    new.names = paste0("A_", Cells(x = pbmc_small[["pca"]]))
head(x = Cells(x = renamed.dimreduc))
# Rename cells in a Seurat object
head(x = colnames(x = pbmc_small))
pbmc_small <- RenameCells(object = pbmc_small, add.cell.id = "A")</pre>
head(x = colnames(x = pbmc_small))
```

RowMergeSparseMatrices

Merge Sparse Matrices by Row

Description

Merge two or more sparse matrices by rowname.

Usage

```
RowMergeSparseMatrices(mat1, mat2)
```

Arguments

mat1 First matrix

mat2 Second matrix or list of matrices

Details

Shared matrix rows (with the same row name) will be merged, and unshared rows (with different names) will be filled with zeros in the matrix not containing the row.

SaveSeuratRds 99

Value

Returns a sparse matrix

SaveSeuratRds

Save and Load Seurat Objects from Rds files

Description

Save and Load Seurat Objects from Rds files

Usage

```
SaveSeuratRds(
  object,
  file = NULL,
  move = TRUE,
  destdir = deprecated(),
  relative = FALSE,
  ...
)
LoadSeuratRds(file, ...)
```

Arguments

```
A Seurat object
object
file
                  Path to save object to; defaults to file.path(getwd(), paste0(Project(object),
                   ".Rds"))
                  Move on-disk layers into dirname(file)
move
destdir
                  [Deprecated]
relative
                  Save relative paths instead of absolute ones
                  Arguments passed on to base::saveRDS, base::readRDS
                  ascii a logical. If TRUE or NA, an ASCII representation is written; otherwise
                       (default), a binary one is used. See the comments in the help for save.
                  version the workspace format version to use. NULL specifies the current default
                       version (3). The only other supported value is 2, the default from R 1.4.0
                       to R 3.5.0.
                  compress a logical specifying whether saving to a named file is to use "gzip"
                       compression, or one of "gzip", "bzip2" or "xz" to indicate the type of
                       compression to be used. Ignored if file is a connection.
                  refhook a hook function for handling reference objects.
```

Value

Invisibly returns file

100 SaveSeuratRds

Progress Updates with progressr

This function uses **progressr** to render status updates and progress bars. To enable progress updates, wrap the function call in with_progress or run handlers(global = TRUE) before running this function. For more details about **progressr**, please read vignette("progressr-intro")

Note

This function requires the fs package to be installed

See Also

```
saveRDS() readRDS()
```

Examples

```
if (requireNamespace("fs", quietly = TRUE)) {
 # Write out with DelayedArray
 if (requireNamespace("HDF5Array", quietly = TRUE)) {
    pbmc <- pbmc_small</pre>
   pbmc[["disk"]] <- CreateAssay50bject(list(</pre>
      mem = LayerData(pbmc, "counts"),
      disk = as(LayerData(pbmc, "counts"), "HDF5Array")
   ))
    # Save `pbmc` to an Rds file
    out <- tempfile(fileext = ".Rds")</pre>
    SaveSeuratRds(pbmc, file = out)
    # Object cache
    obj <- readRDS(out)</pre>
    Tool(obj, "SaveSeuratRds")
    # Load the saved object with on-disk layers back into memory
    pbmc2 <- LoadSeuratRds(out)</pre>
   pbmc2
   pbmc2[["disk"]]
 }
 # Write out with BPCells
 if (requireNamespace("BPCells", quietly = TRUE)) {
    pbmc <- pbmc_small</pre>
    bpm <- BPCells::write_matrix_dir(LayerData(pbmc, "counts"), dir = tempfile())</pre>
   bph <- BPCells::write_matrix_hdf5(</pre>
      LayerData(pbmc, "counts"),
      path = tempfile(fileext = ".h5"),
      group = "counts"
   pbmc[["disk"]] <- CreateAssay50bject(list(dir = bpm, h5 = bph))</pre>
    # Save `pbmc` to an Rds file
```

Segmentation-class 101

```
out <- tempfile(fileext = ".Rds")
SaveSeuratRds(pbmc, file = out)

# Object cache
obj <- readRDS(out)
Tool(obj, "SaveSeuratRds")

# Load the saved object with on-disk layers back into memory
pbmc2 <- LoadSeuratRds(out)
pbmc2
pbmc2[["disk"]]
}</pre>
```

Segmentation-class

The Segmentation Class

Description

The Segmentation Class

See Also

Segmentation methods: Segmentation-methods

 $Segmentation\ layer\ classes:\ Centroids-methods,\ Molecules-class,\ Molecules-methods,\ Segmentation-methods$

Segmentation-methods Segmentation *Methods*

Description

Methods for Segmentation objects

```
## S3 method for class 'Segmentation'
Cells(x, ...)

## S3 method for class 'Segmentation'
GetTissueCoordinates(object, full = TRUE, ...)

## S3 method for class 'Segmentation'
RenameCells(object, new.names = NULL, ...)
```

```
## S3 method for class 'Segmentation'
lengths(x, use.names = TRUE)

## S3 method for class 'Segmentation'
subset(x, cells = NULL, ...)

## S4 method for signature 'Segmentation, ANY, ANY, ANY'
x[i, j, ..., drop = TRUE]

## S4 method for signature 'Segmentation'
coordinates(obj, full = TRUE, ...)

## S4 method for signature 'Segmentation'
show(object)
```

Arguments

x, object, obj A Segmentation object
... Arguments passed to other methods

full Expand the coordinates to the full polygon

new.names vector of new cell names

use.names Ignored

i, cells A vector of cells to keep; if NULL, defaults to all cells

j, drop Ignored

Details

Cells: Get cell names

GetTissueCoordinates, coordinates: Get tissue coordinates

RenameCells: Update cell names

lengths: Generate a run-length encoding of the cells present
subset, [: Subset a Segmentation object to certain cells
show: Display an object summary to stdout

Value

Cells: A vector of cell names

GetTissueCoordinates, coordinates: A data frame with three columns:

- "x": the x-coordinate
- "y": the y-coordinate
- "cell" or "ID": the cell name

set-if-null 103

If full is TRUE, then each coordinate will indicate a vertex for the cell polygon; otherwise, each coordinate will indicate a centroid for the cell. Note: GetTissueCoordinates

RenameCells: object with the cells renamed to new.names

lengths: An rle object for the cells

subset, [: x subsetted to the cells specified by cells/i

show: Invisibly returns NULL

Progress Updates with progressr

The following methods use **progressr** to render status updates and progress bars:

• RenameCells

To enable progress updates, wrap the function call in with_progress or run handlers(global = TRUE) before running this function. For more details about **progressr**, please read **vignette("progressr-intro")**

Parallelization with future

The following methods use **future** to enable parallelization:

• RenameCells

Parallelization strategies can be set using plan. Common plans include "sequential" for non-parallelized processing or "multisession" for parallel evaluation using multiple R sessions; for other plans, see the "Implemented evaluation strategies" section of ?future::plan. For a more thorough introduction to **future**, see vignette("future-1-overview")

See Also

Segmentation-class

 $Segmentation\ layer\ classes:\ Centroids-methods,\ Molecules-class,\ Molecules-methods,\ Segmentation-class$

set-if-null

Set If or If Not NULL

Description

Set a default value depending on if an object is NULL

Usage

```
x %||% y
```

x %iff% y

104 Seurat-class

Arguments

x An object to testy A default value

Value

```
For %||%: y if x is NULL; otherwise x
For %iff%: y if x is not NULL; otherwise x
```

Author(s)

For %||%: rlang developers

See Also

```
rlang::%||%
```

Examples

```
# Set if NULL
1 %||% 2
NULL %||% 2

# Set if *not* NULL
1 %iff% 2
NULL %iff% 2
```

Seurat-class

The Seurat Class

Description

The Seurat object is a representation of single-cell expression data for R; each Seurat object revolves around a set of cells and consists of one or more Assay objects, or individual representations of expression data (eg. RNA-seq, ATAC-seq, etc). These assays can be reduced from their high-dimensional state to a lower-dimension state and stored as DimReduc objects. Seurat objects also store additional metadata, both at the cell and feature level (contained within individual assays). The object was designed to be as self-contained as possible, and easily extendable to new methods.

Slots

```
assays A list of assays for this project
meta.data Contains meta-information about each cell, starting with number of features detected
(nFeature) and the original identity class (orig.ident); more information is added using AddMetaData
active.assay Name of the active, or default, assay; settable using DefaultAssay
active.ident The active cluster identity for this Seurat object; settable using Idents
```

Seurat-validity 105

```
graphs A list of Graph objects
neighbors ...
reductions A list of dimensional reduction objects for this object
images A list of spatial image objects
project.name Name of the project
misc A list of miscellaneous information
version Version of Seurat this object was built under
commands A list of logged commands run on this Seurat object
tools A list of miscellaneous data generated by other tools, should be filled by developers only
using Tool<-
```

See Also

```
Seurat object, validity, and interaction methods $.Seurat(), Seurat-validity, [[.Seurat(),
[[<-,Seurat,NULL,[[<-,Seurat,dim.Seurat(),dimnames.Seurat(),merge.Seurat(),names.Seurat(),
subset.Seurat()</pre>
```

Seurat-validity

Seurat Object Validity

Description

Validation of Seurat objects is handled by validObject

See Also

```
validObject
```

Seurat object, validity, and interaction methods \$.Seurat(), Seurat-class, [[.Seurat(), [[<-, Seurat, NULL, [[<-, Seurat, dim.Seurat(), dimnames.Seurat(), merge.Seurat(), names.Seurat(), subset.Seurat()

SeuratCommand-class

The SeuratCommand Class

Description

The SeuratCommand is used for logging commands that are run on a Seurat object; it stores parameters and timestamps

Slots

```
name Command name
time.stamp Timestamp of when command was tun
assay.used Optional name of assay used to generate SeuratCommand object
call.string String of the command call
params List of parameters used in the command call
```

106 Simplify

See Also

Command log object and interaction methods \$.SeuratCommand(), .DollarNames.SeuratCommand(), LogSeuratCommand(), [.SeuratCommand(), as.list.SeuratCommand()

show, LogMap-method

LogMap Object Overview

Description

Overview of a LogMap object

Usage

```
## S4 method for signature 'LogMap'
show(object)
```

Arguments

object

A LogMap object

Value

Prints summary to stdout and invisibly returns NULL

Simplify

Simplify Geometry

Description

Simplify Geometry

Simplify segmentations by reducing the number of vertices

Usage

```
Simplify(coords, tol, topologyPreserve = TRUE)
## S3 method for class 'Spatial'
Simplify(coords, tol, topologyPreserve = TRUE)
```

Arguments

coords

A 'Segmentation' object

tol

Numerical tolerance value to be used by the Douglas-Peuker algorithm

topologyPreserve

Logical determining if the algorithm should attempt to preserve the topology of

the original geometry

SpatialImage-class 107

Value

A simplified version of coords

A 'Segmentation' object with simplified segmentation vertices

SpatialImage-class The SpatialImage class

Description

The SpatialImage class is a virtual class representing spatial information for Seurat. All spatial image information must inherit from this class for use with Seurat objects

Slots

assay Name of assay to associate image data with; will give this image priority for visualization when the assay is set as the active/default assay in a Seurat object

key A one-length character vector with the object's key; keys must be one or more alphanumeric characters followed by an underscore "_" (regex pattern "^[a-zA-Z][a-zA-Z0-9]*_\$")

See Also

SpatialImage-methods for a list of required and provided methods

SpatialImage-methods SpatialImage methods

Description

Methods defined on the SpatialImage class. Some of these methods must be overridden in order to ensure proper functionality of the derived classes (see **Required methods** below). Other methods are designed to work across all SpatialImage-derived subclasses, and should only be overridden if necessary

```
## S3 method for class 'SpatialImage'
Cells(x, ...)
## S3 method for class 'SpatialImage'
DefaultAssay(object, ...)
## S3 replacement method for class 'SpatialImage'
DefaultAssay(object, ...) <- value</pre>
```

```
## S3 method for class 'SpatialImage'
GetImage(object, mode = c("grob", "raster", "plotly", "raw"), ...)
## S3 method for class 'SpatialImage'
GetTissueCoordinates(object, ...)
## S3 method for class 'SpatialImage'
IsGlobal(object, ...)
## S3 method for class 'SpatialImage'
Key(object, ...)
## S3 replacement method for class 'SpatialImage'
Key(object, ...) <- value</pre>
## S3 method for class 'SpatialImage'
Radius(object)
## S3 method for class 'SpatialImage'
RenameCells(object, new.names = NULL, ...)
## S3 method for class 'SpatialImage'
x[i, ...]
## S3 method for class 'SpatialImage'
dim(x)
## S3 method for class 'SpatialImage'
subset(x, cells, ...)
## S4 method for signature 'SpatialImage'
show(object)
```

Arguments

x, object A SpatialImage-derived object... Arguments passed to other methods

value Depends on the method:

DefaultAssay<- Assay that the image should be associated with

Key<- New key for the image

mode How to return the image; should accept one of "grob", "raster", "plotly", or

"raw"

new.names vector of new cell names
i, cells A vector of cells to keep

Value

[Override] Cells: should return cell names

SpatialImage-methods 109

DefaultAssay: The associated assay of a SpatialImage-derived object

DefaultAssay<-: object with the associated assay updated

[Override] GetImage: The image data from a SpatialImage-derived object

[Override] GetTissueCoordinates: ...

IsGlobal: returns TRUE as images are, by default, global

Key: The key for a SpatialImage-derived object

Key<-: object with the key set to value

Radius: The spot radius size; by default, returns NULL

[Override] RenameCells: object with the new cell names

[, subset: x/object for only the cells requested

[Override] dim: The dimensions of the image data in (Y, X) format

show: Prints summary to stdout and invisibly returns NULL

Functions

- Cells(SpatialImage): Get the cell names from an image ([Override])
- DefaultAssay(SpatialImage): Get the associated assay of a SpatialImage-derived object
- DefaultAssay(SpatialImage) <- value: Set the associated assay of a SpatialImage-derived object
- GetImage(SpatialImage): Get the image data from a SpatialImage-derived object
- GetTissueCoordinates(SpatialImage): Get tissue coordinates for a SpatialImage-derived object ([Override])
- IsGlobal(SpatialImage): Globality test for SpatialImage-derived object
- Key(SpatialImage): Get the key for a SpatialImage-derived object
- Key(SpatialImage) <- value: Set the key for a SpatialImage-derived object
- Radius(SpatialImage): Get the spot radius size
- RenameCells(SpatialImage): Rename cells in a SpatialImage-derived object ([Override])
- [: Subset a SpatialImage-derived object
- dim(SpatialImage): Get the plotting dimensions of an image ([Override])
- subset(SpatialImage): Subset a SpatialImage-derived object ([Override])
- show(SpatialImage): Overview of a SpatialImage-derived object

Provided methods

These methods are defined on the SpatialImage object and should not be overridden without careful thought

- DefaultAssay and DefaultAssay<-
- Key and Key<-
- GetImage; this method *can* be overridden to provide image data, normally returns empty image data. If overridden, should default to returning a grob object

110 split.Assay

- IsGlobal
- Radius; this method can be overridden to provide a spot radius for image objects
- [; this method *can* be overridden to change default subset behavior, normally returns subset(x = x, cells = i). If overridden, should only accept i

Required methods

All subclasses of the SpatialImage class must define the following methods; simply relying on the SpatialImage method will result in errors. For required parameters and their values, see the Usage and Arguments sections

Cells Return the cell/spot barcodes associated with each position

dim Return the dimensions of the image for plotting in (Y, X) format

GetTissueCoordinates Return tissue coordinates; by default, must return a two-column data. frame with x-coordinates in the first column and y-coordinates in the second

Radius Return the spot radius; returns NULL by default for use with non-spot image technologies

RenameCells Rename the cell/spot barcodes for this image

subset Subset the image data by cells/spots

These methods are used throughout Seurat, so defining them and setting the proper defaults will allow subclasses of SpatialImage to work seamlessly

See Also

```
DefaultAssay
GetImage
GetTissueCoordinates
IsGlobal
Key
RenameCells
```

split.Assay

Split an Assay

Description

```
Split an Assay
```

Usage

```
## S3 method for class 'Assay'
split(x, f, drop = FALSE, layers = NA, ...)
```

split.Assay5

Arguments

Χ	An Assay object
f	a 'factor' in the sense that as.factor(f) defines the grouping, or a list of such factors in which case their interaction is used for the grouping. If x is a data frame, f can also be a formula of the form \sim g to split by the variable g, or more generally of the form \sim g1 + + gk to split by the interaction of the variables g1, , gk, where these variables are evaluated in the data frame x using the usual non-standard evaluation rules.
drop	logical indicating if levels that do not occur should be dropped (if f is a factor or a list).
layers	Names of layers to include in the split; pass NA for all layers; pass NULL for the default layer
	Ignored

Value

Returns a v5 assay with splitted layers

See Also

```
v3 Assay object, validity, and interaction methods: $.Assay(), Assay-class, Assay-validity, CreateAssayObject(), [.Assay(), [.Assay(), dim.Assay(), dimnames.Assay(), merge.Assay(), subset.Assay()
```

split.Assay5

Split an Assay

Description

Split an Assay

Usage

```
## S3 method for class 'Assay5'
split(
    x,
    f,
    drop = FALSE,
    layers = c("counts", "data"),
    ret = c("assay", "multiassays", "layers"),
    ...
)
```

112 split.Assay5

Arguments

x	An Assay5 object
f	a 'factor' in the sense that <code>as.factor(f)</code> defines the grouping, or a list of such factors in which case their interaction is used for the grouping. If x is a data frame, f can also be a formula of the form $\sim g$ to split by the variable g , or more generally of the form $\sim g1 + \ldots + gk$ to split by the interaction of the variables $g1, \ldots, gk$, where these variables are evaluated in the data frame x using the usual non-standard evaluation rules.
drop	logical indicating if levels that do not occur should be dropped (if f is a factor or a list).
layers	Names of layers to include in the split; pass NA for all layers; pass NULL for the default layer
ret	Type of return value; choose from:
	• "assay": a single Assay5 object
	• "multiassay": a list of Assay5 objects
	"layers": a list of layer matrices
	Ignored

Value

Depends on the value of ret:

- "assay": x with the layers requested in layers split based on f; all other layers are left as-is
- "multiassay": a list of Assay5 objects; the list contains one value per split and each assay contains only the layers requested in layers with the key set to the split
- "layers": a list of matrices of length length(assays) * length(unique(f)); the list is named as "layer.split"

Progress Updates with progressr

This function uses **progressr** to render status updates and progress bars. To enable progress updates, wrap the function call in with_progress or run handlers(global = TRUE) before running this function. For more details about **progressr**, please read vignette("progressr-intro")

See Also

```
v5 Assay object, validity, and interaction methods: $.Assay5(), Assay5-class, Assay5-validity, [.Assay5(), [[.Assay5(), dim.Assay5(), dimnames.Assay5(), merge.Assay5(), subset.Assay5()
```

Stdev 113

Stdev

Get the standard deviations for an object

Description

Get the standard deviations for an object

Usage

```
Stdev(object, ...)
## S3 method for class 'DimReduc'
Stdev(object, ...)
## S3 method for class 'Seurat'
Stdev(object, reduction = "pca", ...)
```

Arguments

object An object

... Arguments passed to other methods

reduction Name of reduction to use

Value

The standard deviations

Examples

```
# Get the standard deviations for each PC from the DimReduc object
Stdev(object = pbmc_small[["pca"]])
# Get the standard deviations for each PC from the Seurat object
Stdev(object = pbmc_small, reduction = "pca")
```

StitchMatrix

Stitch Matrices Together

Description

Stitch Matrices Together

Usage

```
StitchMatrix(x, y, rowmap, colmap, ...)
```

114 subset.Assay

Arguments

x A matrix

y One or more matrices of the same class or coercible to the same class as x

rowmap, colmap LogMaps describing the row and cell membership of each matrix; the LogMap entries are assumed to be in the order of c(x, y)

... Arguments passed to other methods

Value

A single matrix of type class(x) consisting of all values in component matrices

subset. Assay Subset an Assay

Description

Subset an Assay

Usage

```
## S3 method for class 'Assay'
subset(x, cells = NULL, features = NULL, ...)
```

Arguments

```
x An Assay objectcells Cell namesfeatures Feature names... Ignored
```

Value

x with just the cells and features specified by cells and features

See Also

```
v3 Assay object, validity, and interaction methods: $.Assay(), Assay-class, Assay-validity, CreateAssayObject(), [.Assay(), [.Assay(), dim.Assay(), dimnames.Assay(), merge.Assay(), split.Assay()
```

```
rna <- pbmc_small[["RNA"]]
rna2 <- subset(rna, features = VariableFeatures(rna))
rna2</pre>
```

subset.Assay5

subset.Assay5

Subset an Assay

Description

Subset an Assay

Usage

```
## S3 method for class 'Assay5'
subset(x, cells = NULL, features = NULL, layers = NULL, ...)
```

Arguments

x An Assay5 object

cells Cell names
features Feature names

layers Layer to keep; defaults to all layers

... Ignored

Value

x with just the cells and features specified by cells and features for the layers specified by layers

See Also

```
v5 Assay object, validity, and interaction methods: $.Assay5(), Assay5-class, Assay5-validity, [.Assay5(), [[.Assay5(), dim.Assay5(), dimnames.Assay5(), merge.Assay5(), split.Assay5()
```

subset.DimReduc

Subset a Dimensional Reduction

Description

Subset a DimReduc object

Usage

```
## S3 method for class 'DimReduc'
subset(x, cells = NULL, features = NULL, ...)
```

116 subset.Seurat

Arguments

```
x A DimReduc object cells, features

Cells and features to keep during the subset

Ignored
```

Value

x for cells cells and features features

See Also

Dimensional reduction object, validity, and interaction methods CreateDimReducObject(), DimReduc-class, DimReduc-validity, [.DimReduc(), [[.DimReduc(), dim.DimReduc(), merge.DimReduc(), print.DimReduc()

subset.Seurat

Subset Seurat Objects

Description

Subset Seurat Objects

Usage

```
## S3 method for class 'Seurat'
subset(
    x,
    subset,
    cells = NULL,
    features = NULL,
    idents = NULL,
    return.null = FALSE,
    ...
)

## S3 method for class 'Seurat'
x[i, j, ...]
```

Arguments

```
x A Seurat object
subset Logical expression indicating features/variables to keep
cells, j A vector of cell names or indices to keep
features, i A vector of feature names or indices to keep
idents A vector of identity classes to keep
return.null If no cells are requested, return a NULL; by default, throws an error
... Arguments passed to WhichCells
```

Theta 117

Value

```
subset: A subsetted Seurat object
[: object x with features i and cells j
```

See Also

```
WhichCells
```

```
Seurat object, validity, and interaction methods $.Seurat(), Seurat-class, Seurat-validity,
[[.Seurat(), [[<-,Seurat, NULL, [[<-,Seurat, dim.Seurat(), dimnames.Seurat(), merge.Seurat(),
names.Seurat()</pre>
```

Examples

```
# `subset` examples
subset(pbmc_small, subset = MS4A1 > 4)
subset(pbmc_small, subset = `DLGAP1-AS1` > 2)
subset(pbmc_small, idents = '0', invert = TRUE)
subset(pbmc_small, subset = MS4A1 > 3, slot = 'counts')
subset(pbmc_small, features = VariableFeatures(object = pbmc_small))
# `[` examples
pbmc_small[VariableFeatures(object = pbmc_small), ]
pbmc_small[, 1:10]
```

Theta

Get the offset angle

Description

Get the offset angle

Usage

Theta(object)

Arguments

object An object

Tool

Tool

Get and Set Additional Tool Data

Description

Use Tool to get tool data. If no additional arguments are provided, will return a vector with the names of tools in the object.

Usage

```
Tool(object, ...)
Tool(object, ...) <- value
## S3 method for class 'Seurat'
Tool(object, slot = NULL, ...)
## S3 replacement method for class 'Seurat'
Tool(object, ...) <- value</pre>
```

Arguments

object	An object
	Arguments passed to other methods
value	Information to be added to tool list
slot	Name of tool to pull

Value

If no additional arguments, returns the names of the tools in the object; otherwise returns the data placed by the tool requested

Note

For developers: set tool data using Tool<-. Tool<- will automatically set the name of the tool to the function that called Tool<-, so each function gets one entry in the tools list and cannot overwrite another function's entry. The automatic naming will also remove any method identifiers (eg. RunPCA. Seurat will become RunPCA); please plan accordingly

```
# Example function that adds unstructured data to tools
MyTool <- function(object) {
   sample.tool.output <- matrix(rnorm(n = 16), nrow = 4)
   # Note: `Tool<-` must be called from within a function
   # and the name of the tool will be generated from the function name
   Tool(object) <- sample.tool.output</pre>
```

UpdateSeuratObject 119

```
return(object)
}

# Run our tool
set.seed(42L)
pbmc_small <- MyTool(pbmc_small)

# Get a list of tools run
Tool(pbmc_small)

# Access specific tool data
Tool(pbmc_small, slot = "MyTool")</pre>
```

UpdateSeuratObject

Update old Seurat object to accommodate new features

Description

Updates Seurat objects to new structure for storing data/calculations. For Seurat v3 objects, will validate object structure ensuring all keys and feature names are formed properly.

Usage

```
UpdateSeuratObject(object)
```

Arguments

object

Seurat object

Value

Returns a Seurat object compatible with latest changes

```
## Not run:
updated_seurat_object = UpdateSeuratObject(object = old_seurat_object)
## End(Not run)
```

120 Version

UpdateSlots

Update slots in an object

Description

Update slots in an object

Usage

```
UpdateSlots(object)
```

Arguments

object

An object to update

Value

object with the latest slot definitions

Version

Get Version Information

Description

Get Version Information

Usage

```
Version(object, ...)
## S3 method for class 'Seurat'
Version(object, ...)
```

Arguments

object An object

... Arguments passed to other methods

```
Version(pbmc_small)
```

WhichCells 121

WhichCells

Identify cells matching certain criteria

Description

Returns a list of cells that match a particular set of criteria such as identity class, high/low values for particular PCs, etc.

Usage

```
WhichCells(object, ...)
## S3 method for class 'Assay'
WhichCells(object, cells = NULL, expression, invert = FALSE, ...)
## S3 method for class 'Seurat'
WhichCells(
  object,
  cells = NULL,
  idents = NULL,
  expression,
  slot = "data",
  invert = FALSE,
  downsample = Inf,
  seed = 1,
  ...
)
```

Arguments

object	An object
	Arguments passed on to CellsByIdentities
	return.null If no cells are requested, return a NULL; by default, throws an error
cells	Subset of cell names
expression	A predicate expression for feature/variable expression, can evaluate anything that can be pulled by FetchData; please note, you may need to wrap feature names in backticks (``) if dashes between numbers are present in the feature name
invert	Invert the selection of cells
idents	A vector of identity classes to keep
slot	Slot to pull feature data for
downsample	Maximum number of cells per identity class, default is Inf; downsampling will happen after all other operations, including inverting the cell selection
seed	Random seed for downsampling. If NULL, does not set a seed

[.Assay

Value

A vector of cell names

See Also

FetchData

Examples

```
WhichCells(pbmc_small, idents = 2)
WhichCells(pbmc_small, expression = MS4A1 > 3)
levels(pbmc_small)
WhichCells(pbmc_small, idents = c(1, 2), invert = TRUE)
```

[.Assay

Layer Data

Description

Get and set layer data

Usage

```
## S3 method for class 'Assay'
x[i = missing_arg(), j = missing_arg(), ...]
## S4 replacement method for signature 'Assay, character, ANY, ANY'
x[i, j, ...] <- value</pre>
```

Arguments

X	An Assay object
i	Name of layer data to get or set
j	Ignored
	Arguments passed to LayerData
value	A matrix-like object to add as a new layer

Value

```
[: The layer data for layer i
[<-: x with layer data value saved as i
```

[.Assay5

See Also

```
LayerData
```

```
v3 Assay object, validity, and interaction methods: $.Assay(), Assay-class, Assay-validity, CreateAssayObject(), [[.Assay(), dim.Assay(), dimnames.Assay(), merge.Assay(), split.Assay(), subset.Assay()
```

Examples

```
rna <- pbmc_small[["RNA"]]
# Get a vector of layer names in this assay
rna[]
# Fetch layer data
rna["data"][1:10, 1:4]
# Set layer data
rna["data"] <- rna["counts"]
rna["data"][1:10, 1:4]</pre>
```

[.Assay5

Layer Data

Description

Get and set layer data

Usage

```
## S3 method for class 'Assay5'
x[i = missing_arg(), j = missing_arg(), ...]
## S4 replacement method for signature 'Assay5, character, ANY, ANY'
x[i, j, ...] <- value</pre>
```

Arguments

X	An Assay5 object
i	Name of layer data to get or set
j	Ignored
	Arguments passed to LayerData
value	A matrix-like object to add as a new layer

124 [.DimReduc

Value

```
[: The layer data for layer i
[<-: x with layer data value saved as i
```

See Also

LayerData

```
v5 Assay object, validity, and interaction methods: $.Assay5(), Assay5-class, Assay5-validity, [[.Assay5(), dim.Assay5(), dimnames.Assay5(), merge.Assay5(), split.Assay5(), subset.Assay5()
```

[.DimReduc

Get Feature Loadings

Description

Pull feature loadings from a dimensional reduction

Usage

```
## S3 method for class 'DimReduc'
x[i, j, drop = FALSE, ...]
```

Arguments

Х	A DimReduc object
i	Feature identifiers or indices
j	Dimension identifiers or indices
drop	Coerce the result to the lowest possible dimension; see drop for further details
	Arguments passed to other methods

Details

[does not distinguish between projected and unprojected feature loadings; to select whether projected or unprojected loadings should be pulled, please use Loadings

Value

Feature loadings for features i and dimensions j

See Also

Loadings

```
Dimensional reduction object, validity, and interaction methods CreateDimReducObject(), DimReduc-class, DimReduc-validity, [[.DimReduc(), dim.DimReduc(), merge.DimReduc(), print.DimReduc(), subset.DimReduc()
```

[.SeuratCommand 125

Examples

```
pca <- pbmc_small[["pca"]]
pca[1:10, 1:5]</pre>
```

[.SeuratCommand

Command Log Data Access

Description

Access data from a SeuratCommand object

Usage

```
## S3 method for class 'SeuratCommand' x[i, \ldots]
```

Arguments

x A SeuratCommand object

i The name of a command log slot

... Ignored

Value

[: Slot i from x

See Also

Command log object and interaction methods \$. SeuratCommand(), .DollarNames.SeuratCommand(), LogSeuratCommand(), SeuratCommand-class, as.list.SeuratCommand()

```
cmd <- pbmc_small[["NormalizeData.RNA"]]
cmd["call.string"]</pre>
```

[[.Assay

[[.Assay

Feature-Level Meta Data

Description

Get and set feature-level meta data

Usage

```
## S3 method for class 'Assay'
x[[i, ..., drop = FALSE]]

## S4 replacement method for signature 'Assay,ANY,ANY,ANY'
x[[i, j, ...]] <- value

## S3 method for class 'Assay'
head(x, n = 10L, ...)

## S3 method for class 'Assay'
tail(x, n = 10L, ...)

## S4 replacement method for signature 'Assay,missing,missing,data.frame'
x[[i, j, ...]] <- value</pre>
```

Arguments

X	An Assay object
i	Name of feature-level meta data to fetch or add
	Ignored
drop	See drop
j	Ignored
value	Feature-level meta data to add
n	Number of meta data rows to show

Value

```
[[: The feature-level meta data for i
[[<-: x with value added as i in feature-level meta data
head: The first n rows of feature-level meta data
tail: the last n rows of feature-level meta data
```

See Also

```
v3 Assay object, validity, and interaction methods: $.Assay(), Assay-class, Assay-validity, CreateAssayObject(), [.Assay(), dim.Assay(), dimnames.Assay(), merge.Assay(), split.Assay(), subset.Assay()
```

[[.Assay5

Examples

```
rna <- pbmc_small[["RNA"]]

# Pull the entire feature-level meta data frame
head(rna[[]])

# Pull a specific column of feature-level meta data
head(rna[["vst.mean"]])
head(rna[["vst.mean", drop = TRUE]])

# `head` and `tail` can be used to quickly view feature-level meta data
head(rna)

tail(rna)</pre>
```

[[.Assay5

Feature-Level Meta Data

Description

Get and set feature-level meta data

Usage

```
## S3 method for class 'Assay5'
x[[i, j, ..., drop = FALSE]]

## S4 replacement method for signature 'Assay5,ANY,ANY,ANY'
x[[i, j, ...]] <- value

## S3 method for class 'Assay5'
head(x, n = 10L, ...)

## S3 method for class 'Assay5'
tail(x, n = 10L, ...)</pre>
```

Arguments

X	An Assay5 object
i	Name of feature-level meta data to fetch or add
j	Ignored
• • •	Ignored
drop	See drop
value	Feature-level meta data to add
n	Number of meta data rows to show

128 [[.DimReduc

Value

```
[[: The feature-level meta data for i
```

[[<-: x with value added as i in feature-level meta data

head: The first n rows of feature-level meta data

tail: the last n rows of feature-level meta data

See Also

```
v5 Assay object, validity, and interaction methods: $.Assay5(), Assay5-class, Assay5-validity, [.Assay5(), dim.Assay5(), dimnames.Assay5(), merge.Assay5(), split.Assay5(), subset.Assay5()
```

[[.DimReduc

Get Cell Embeddings

Description

Pull cell embeddings from a dimensional reduction

Usage

```
## S3 method for class 'DimReduc'
x[[i, j, drop = FALSE, ...]]
```

Arguments

X	A DimReduc object
i	Cell names or indices

j Dimension identifiers or indices

drop Coerce the result to the lowest possible dimension; see drop for further details

... Arguments passed to other methods

Value

Cell embeddings for cells i and dimensions j

See Also

Embeddings

```
Dimensional reduction object, validity, and interaction methods CreateDimReducObject(), DimReduc-class, DimReduc-validity, [.DimReduc(), dim.DimReduc(), merge.DimReduc(), print.DimReduc(), subset.DimReduc()
```

```
pca <- pbmc_small[["pca"]]
pca[[1:10, 1:5]]</pre>
```

[[.Seurat 129

[[.Seurat

Subobjects and Cell-Level Meta Data

Description

The [[operator pulls either subobjects (eg. v3 or v5 assays, dimensional reduction information, or nearest-neighbor graphs) or cell-level meta data from a Seurat object

Usage

```
## $3 method for class 'Seurat'
x[[i = missing_arg(), ..., drop = FALSE, na.rm = FALSE]]
## $3 method for class 'Seurat'
head(x, n = 10L, ...)
## $3 method for class 'Seurat'
tail(x, n = 10L, ...)
```

Arguments

X	A Seurat object
i	Name of cell-level meta data
	Ignored
drop	See drop
na.rm	Remove cells where meta data is all NA
n	Number of meta data rows to show

Value

Varies based on the value of i:

- If i is missing, a data frame with cell-level meta data
- If i is a vector with cell-level meta data names, a data frame (or vector of drop = TRUE) with cell-level meta data requested
- If i is a one-length character with the name of a subobject, the subobject specified by i

head: The first n rows of cell-level metadata tail: The last n rows of cell-level metadata

See Also

See here for adding meta data with [[<-, here for adding subobjects with [[<-, and here for removing subobjects and cell-level meta data with [[<-

```
Seurat object, validity, and interaction methods $.Seurat(), Seurat-class, Seurat-validity,
[[<-,Seurat,NULL,[[<-,Seurat,dim.Seurat(),dimnames.Seurat(),merge.Seurat(),names.Seurat(),
subset.Seurat()</pre>
```

[[<-,Seurat

Examples

```
# Get the cell-level metadata data frame
head(pbmc_small[[]])

# Pull specific metadata information
head(pbmc_small[[c("letter.idents", "groups")]])
head(pbmc_small[["groups", drop = TRUE]])

# Get a sub-object (eg. an `Assay` or `DimReduc`)
pbmc_small[["RNA"]]
pbmc_small[["pca"]]

# Get the first 10 rows of cell-level metadata
head(pbmc_small)

# Get the last 10 rows of cell-level metadata
tail(pbmc_small)
```

[[<-,Seurat</pre>

Add Subobjects

Description

Add subobjects containing expression, dimensional reduction, or other containerized data to a Seurat object. Subobjects can be accessed with [[and manipulated directly within the Seurat object or used independently

Usage

```
## S4 replacement method for signature 'Seurat,character,missing,Assay'
x[[i, j, ...]] <- value

## S4 replacement method for signature 'Seurat,character,missing,Assay5'
x[[i, j, ...]] <- value

## S4 replacement method for signature 'Seurat,character,missing,DimReduc'
x[[i, j, ...]] <- value

## S4 replacement method for signature 'Seurat,character,missing,Graph'
x[[i, j, ...]] <- value

## S4 replacement method for signature 'Seurat,character,missing,Neighbor'
x[[i, j, ...]] <- value

## S4 replacement method for signature 'Seurat,character,missing,Neighbor'
x[[i, j, ...]] <- value</pre>
```

[[<-,Seurat,NULL 131

```
## S4 replacement method for signature 'Seurat, character, missing, Spatial Image' x[[i, j, ...]] \leftarrow value
```

Arguments

```
x A Seurat object
i Name to add subobject as
j Ignored
... Ignored
```

value A valid subobject (eg. a v3 or v5 assay, or a dimensional reduction)

Value

x with value added as i

See Also

See here for pulling subobjects using [[, here for adding metadata with [[<-, and here for removing subobjects and cell-level meta data with [[<-

```
Seurat object, validity, and interaction methods $.Seurat(), Seurat-class, Seurat-validity,
[[.Seurat(), [[<-, Seurat, NULL, dim.Seurat(), dimnames.Seurat(), merge.Seurat(), names.Seurat(),
subset.Seurat()</pre>
```

[[<-,Seurat,NULL R

Remove Subobjects and Cell-Level Meta Data

Description

Remove Subobjects and Cell-Level Meta Data

Usage

```
## S4 replacement method for signature 'Seurat, character, missing, `NULL`' x[[i, j, ...]] \leftarrow value
```

Arguments

х	A Seurat object
i	Name(s) of $subobject(s)$ or cell-level meta data to remove
j	Ignored
• • •	Ignored
value	NULL

Value

x with i removed from the object

\$.Assay

See Also

See here for pulling subobjects using [[, here for adding metadata with [[<-, and here for adding subobjects with [[<-

Seurat object, validity, and interaction methods \$.Seurat(), Seurat-class, Seurat-validity,
[[.Seurat(), [[<-, Seurat, dim.Seurat(), dimnames.Seurat(), merge.Seurat(), names.Seurat(),
subset.Seurat()</pre>

\$.Assay

Layer Data

Description

Get and set layer data

Usage

```
## S3 method for class 'Assay'
x$i
## S3 replacement method for class 'Assay'
x$i <- value</pre>
```

Arguments

x An Assay object

i Name of layer data to get or set

value A matrix-like object to add as a new layer

Value

```
$: Layer data for layer i
```

\$<-: x with layer data value saved as i</pre>

See Also

```
v3 Assay object, validity, and interaction methods: Assay-class, Assay-validity, CreateAssayObject(), [.Assay(), [[.Assay(), dim.Assay(), dimnames.Assay(), merge.Assay(), split.Assay(), subset.Assay()
```

\$.Assay5

Examples

```
rna <- pbmc_small[["RNA"]]
# Fetch a layer with `$`
rna$data[1:10, 1:4]
# Add a layer with `$`
rna$data <- rna$counts
rna$data[1:10, 1:4]</pre>
```

\$.Assay5

Layer Data

Description

Get and set layer data

Usage

```
## S3 method for class 'Assay5'
x$i
## S3 replacement method for class 'Assay5'
x$i <- value</pre>
```

Arguments

x An Assay5 object

i Name of layer data to get or set

value A matrix-like object to add as a new layer

Value

```
$: Layer data for layer i
```

\$<-: x with layer data value saved as i</pre>

See Also

```
v5 Assay object, validity, and interaction methods: Assay5-class, Assay5-validity, [.Assay5(), [[.Assay5(), dim.Assay5(), dimnames.Assay5(), merge.Assay5(), split.Assay5(), subset.Assay5()
```

\$.Seurat

\$.Seurat

Cell-Level Meta Data

Description

Get and set cell-level meta data

Usage

```
## S3 method for class 'Seurat'
x$i

## S3 replacement method for class 'Seurat'
x$i, ... <- value

## S4 replacement method for signature 'Seurat, character, missing, data. frame'
x[[i, j, ...]] <- value

## S4 replacement method for signature 'Seurat, missing, missing, data. frame'
x[[i, j, ...]] <- value

## S4 replacement method for signature 'Seurat, character, missing, factor'
x[[i, j, ...]] <- value

## S4 replacement method for signature 'Seurat, character, missing, list'
x[[i, j, ...]] <- value

## S4 replacement method for signature 'Seurat, missing, missing, list'
x[[i, j, ...]] <- value

## S4 replacement method for signature 'Seurat, character, missing, vector'
x[[i, j, ...]] <- value</pre>
```

Arguments

X	A Seurat object
i	Name of cell-level meta data
	Ignored
value	A vector to add as cell-level meta data
j	Ignored

Value

\$: Metadata column i for object x; **note**: unlike [[, \$ drops the shape of the metadata to return a vector instead of a data frame

```
$<-: x with metadata value saved as i</pre>
```

\$.SeuratCommand 135

See Also

Examples

```
# Get metadata using `$'
head(pbmc_small$groups)

# Add metadata using the `$' operator
set.seed(42)
pbmc_small$value <- sample(1:3, size = ncol(pbmc_small), replace = TRUE)
head(pbmc_small[["value"]])</pre>
```

\$.SeuratCommand

Command Log Parameter Access

Description

Pull parameter values from a SeuratCommand object

Usage

```
## S3 method for class 'SeuratCommand' x$i
```

Arguments

- x A SeuratCommand object
- i A parameter name

Value

The value for parameter i

See Also

```
Command log object and interaction methods .DollarNames.SeuratCommand(), LogSeuratCommand(), SeuratCommand-class, [.SeuratCommand(), as.list.SeuratCommand()
```

```
cmd <- pbmc_small[["NormalizeData.RNA"]]
cmd$normalization.method</pre>
```

Index

* assay5	CellsByImage, 23
[.Assay5, 123	Command, 28
[[.Assay5, 127	DefaultAssay, 38
\$.Assay5, 133	Distances, 49
Assay5-class, 15	Embeddings, 51
Assay5-validity, 16	FetchData, 51
CastAssay, 21	GetImage, 57
DefaultLayer, 41	GetTissueCoordinates, 58
dim.Assay5,42	HVFInfo, 59
dimnames.Assay5,46	Images, 65
JoinLayers, 71	Index, 66
merge.Assay5,83	Indices, 66
split.Assay5, 111	IsGlobal, 68
subset.Assay5, 115	Key, 73
* assay	LayerData, 75
[.Assay, 122	Loadings, 77
[[.Assay, 126	Misc, 86
\$.Assay, 132	Stdev, 113
Assay-class, 13	Tool, 118
Assay-validity, 14	Version, 120
CreateAssay5Object, 28	WhichCells, 121
CreateAssayObject, 29	* datasets
dim.Assay,41	pbmc_small, 92
dimnames.Assay,45	* dimnames
merge.Assay,82	Cells, <u>22</u>
split.Assay, 110	dimnames.Assay,45
subset.Assay, 114	dimnames.Assay5,46
* command	dimnames.Seurat,46
$. {\tt DollarNames.SeuratCommand,7}$	* dimreduc
[.SeuratCommand, 125	[.DimReduc, 124
\$.SeuratCommand, 135	[[.DimReduc, 128
as.list.SeuratCommand, 10	CreateDimReducObject, 31
LogSeuratCommand, 81	dim.DimReduc, 43
SeuratCommand-class, 105	DimReduc-class, 47
* data-access	${\tt DimReduc-validity}, 48$
AssayData, 16	merge.DimReduc,84
Assays, 18	print.DimReduc,93
Cells, 22	subset.DimReduc, 115
CellsByIdentities, 23	* fov

FOV-class, 53	UpdateSeuratObject, 119
FOV-methods, 54	* spatialimage
FOV-validity, 57	Radius, 95
* future	SpatialImage-methods, 107
Segmentation-methods, 101	* spatial
* graph	as.Centroids,8
as.Graph, 9	Boundaries, 20
Graph-class, 59	CreateCentroids, 30
* jackstraw	CreateFOV, 32
JackStrawData-methods, 70	CreateMolecules, 33
JS, 72	CreateSegmentation, 34
* logmap	Crop, 37
as.matrix.LogMap, 11	DefaultFOV, 40
droplevels.LogMap, 50	Overlay, 91
intersect.LogMap, 67	Simplify, 106
labels.LogMap, 74	Theta, 117
LogMap, 78	* utils
LogMap-validity, 80	as.sparse, 13
show, LogMap-method, 106	AttachDeps, 19
* neighbor	CheckGC, 27
as.Neighbor, 11	CheckLayersName, 27
Neighbor-methods, 90	DefaultDimReduc, 39
* segmentation	FilterObjects, 53
Centroids-class, 24	IsMatrixEmpty, 68
Centroids-methods, 25	IsNamedList, 69
Molecules-class, 87	PackageCheck, 92
Molecules-methods, 88	RandomName, 95
Segmentation-class, 101	RowMergeSparseMatrices, 98
Segmentation-methods, 101	SaveSeuratRds, 99
* seurat	set-if-null, 103
[[.Seurat, 129	StitchMatrix, 113
[[<-,Seurat, 130	UpdateSlots, 120
[[<-,Seurat,NULL, 131	.DollarNames.JackStrawData
\$.Seurat, 134	(JackStrawData-methods), 70
AddMetaData, 7	.DollarNames.SeuratCommand, 7, 10, 81,
as.Seurat, 12	106, 125, 135
CreateSeuratObject, 35	.FilterObjects, 53
dim.Seurat,44	.onAttach, 19
dimnames.Seurat,46	?future::plan, <i>103</i>
Idents, 62	[, 110
merge.Seurat,84	[,Centroids,character,ANY,ANY-method
names.Seurat, 89	(Centroids-methods), 25
Project, 94	[,Centroids,numeric,ANY,ANY-method
RenameAssays, 96	(Centroids-methods), 25
RenameCells, 97	[,Segmentation,ANY,ANY,ANY-method
Seurat-class, 104	(Segmentation-methods), 101
Seurat-validity, 105	[.Assay, 14, 15, 30, 42, 45, 82, 111, 114, 122
subset.Seurat, 116	126, 132

[.Assay5, 16, 42, 46, 83, 112, 115, 123, 128,	(LogMap), 78
133	<pre>[[<-,LogMap,character,missing,integer-method</pre>
[.DimReduc, 31, 43, 48, 49, 84, 94, 116, 124,	(LogMap), 78
128	<pre>[[<-,LogMap,character,missing,numeric-method</pre>
[.FOV (FOV-methods), 54	(LogMap), 78
[.Seurat (subset.Seurat), 116	<pre>[[<-,Seurat,character,missing,Assay-method</pre>
[.SeuratCommand, 7, 10, 81, 106, 125, 135	([[<-,Seurat), 130]
[.SpatialImage(SpatialImage-methods),	<pre>[[<-,Seurat,character,missing,Assay5-method</pre>
107	([[<-,Seurat), 130]
<pre>[<-,Assay,character,ANY,ANY-method</pre>	<pre>[[<-,Seurat,character,missing,DimReduc-method</pre>
([.Assay), 122	([[<-,Seurat), 130]
<pre>[<-,Assay5,character,ANY,ANY-method</pre>	<pre>[[<-,Seurat,character,missing,Graph-method</pre>
([.Assay5), 123	([[<-, Seurat), 130]
[[, 130	<pre>[[<-,Seurat,character,missing,NULL-method</pre>
[[,LogMap,NULL,missing-method(LogMap),	([[<-,Seurat,NULL),131
78	<pre>[[<-,Seurat,character,missing,Neighbor-method</pre>
[[,LogMap,character,missing-method	([[<-, Seurat), 130]
(LogMap), 78	$\verb [[<-,Seurat,character,missing,SeuratCommand-method \\$
[[,LogMap,missing,missing-method	([[<-, Seurat), 130]
(LogMap), 78	[[<-,Seurat,character,missing,SpatialImage-method
[[.Assay, 14, 15, 30, 42, 45, 82, 111, 114,	([[<-, Seurat), 130
123, 126, 132	[[<-,Seurat,character,missing,data.frame-method
[[.Assay5, 16, 42, 46, 83, 112, 115, 124, 127,	(\$.Seurat), 134
133	[[<-,Seurat,character,missing,factor-method
[[.DimReduc, 31, 43, 48, 49, 84, 94, 116, 124,	(\$.Seurat), 134
128	[[<-,Seurat,character,missing,list-method
[[.FOV (FOV-methods), 54	(\$.Seurat), 134
[[.Seurat, 44, 47, 85, 89, 105, 117, 129, 131, 132, 135]	<pre>[[<-,Seurat,character,missing,vector-method</pre>
[[<-,Seurat, 130	<pre>[[<-,Seurat,missing,missing,data.frame-method</pre>
[[<-,Seurat,NULL, 131	(\$.Seurat), 134
[[<-,Assay,ANY,ANY,ANY-method	<pre>[[<-,Seurat,missing,missing,list-method</pre>
([[.Assay), 126	(\$.Seurat), 134
[[<-,Assay,missing,missing,data.frame-method	[[<seurat ([[<-,seurat),="" 130<="" td=""></seurat>
([[.Assay), 126	\$. Assay, 14, 15, 30, 42, 45, 82, 111, 114, 123,
[[<-,Assay5,ANY,ANY,ANY-method	126, 132
([[.Assay5), 127	\$. Assay5, 16, 42, 46, 83, 112, 115, 124, 128,
[[<-,FOV,character,missing,Centroids-method	133
(FOV-methods), 54	\$.FOV (FOV-methods), 54
<pre>[[<-,FOV,character,missing,Molecules-method</pre>	<pre>\$.JackStrawData</pre>
(FOV-methods), 54	(JackStrawData-methods), 70
[[<-,FOV,character,missing,NULL-method	\$. Seurat, 44, 47, 85, 89, 105, 117, 129, 131,
(FOV-methods), 54	<i>132</i> , 134
[[<-,FOV,character,missing,Segmentation-methology]	
(FOV-methods), 54	\$ <assay (\$.assay),="" 132<="" td=""></assay>
[[<-,LogMap,character,missing,NULL-method	\$< Assay5 (\$. Assay5), 133
(LogMap), 78	\$ <seurat (\$.seurat),="" 134<="" td=""></seurat>
[[<- LogMan character missing character-methol	0%1tt%/set-it-null) 103

AddMataData 7 104	accordinates Commentation mathed
AddMetaData, 7, 104	coordinates, Segmentation-method
AddSamples (merge. Seurat), 84	(Segmentation-methods), 101
as.Centroids, 8	CreateAssay50bject, 28
as.factor, 111, 112	CreateAssayObject, 14, 15, 29, 42, 45, 82,
as. Graph, 9, 59	111, 114, 123, 126, 132
as.list.SeuratCommand, 7, 10, 81, 106, 125,	CreateCentroids, 30
135	CreateDimReducObject, 31, 43, 48, 49, 84,
as.logical.JackStrawData	94, 116, 124, 128
(JackStrawData-methods), 70	CreateFOV, 32
as.matrix.LogMap, 11, 50, 67, 75, 80, 81	CreateMolecules, 33
as.Neighbor, 11	CreateSegmentation, 34
as. Segmentation (as. Centroids), 8	CreateSeuratObject, 35
as. Seurat, 12	Crop, 37
as.sparse, 13, 29	data Carama 15
Assay, 16, 18, 30, 36, 42, 45, 59, 82, 104, 111,	data frame, 15
114, 122, 126, 132	data.frame, 33, 36
Assay (Assay-class), 13	default assay, 17, 47
Assay-class, 13	default layer, 15, 111, 112
Assay-validity, 14	DefaultAssay, 38, 104, 109, 110
Assay5, 28, 29, 42, 46, 83, 112, 115, 123, 127,	DefaultAssay.SpatialImage
133	(SpatialImage-methods), 107
Assay5 (Assay5-class), 15	DefaultAssay<- (DefaultAssay), 38
Assay5-class, 15	DefaultAssay <spatialimage< td=""></spatialimage<>
Assay5-validity, 16	(SpatialImage-methods), 107
AssayData, 16	DefaultBoundary (Boundaries), 20
Assays, 18	DefaultBoundary<- (Boundaries), 20
AttachDeps, 19	DefaultDimReduc, 39
	DefaultFOV, 40
base::readRDS, 99	DefaultFOV<- (DefaultFOV), 40
base::saveRDS, 99	DefaultLayer, 41
Boundaries, 20	DefaultLayer<- (DefaultLayer), 41
	dgCMatrix, 59
CastAssay, 21	dim, <i>110</i>
cat, <i>94</i>	dim. Assay, 14, 15, 30, 41, 45, 82, 111, 114,
Cells, 22, 45–47, 110	123, 126, 132
Cells.Centroids (Centroids-methods), 25	dim. Assay5, 16, 42, 46, 83, 112, 115, 124,
Cells.FOV (FOV-methods), 54	128, 133
Cells.Segmentation	dim.DimReduc, 31, 43, 48, 49, 84, 94, 116,
(Segmentation-methods), 101	124, 128
Cells.SpatialImage	dim.Neighbor(Neighbor-methods),90
(SpatialImage-methods), 107	dim. Seurat, 44, 47, 85, 89, 105, 117, 129,
CellsByIdentities, 23, 121	131, 132, 135
CellsByImage, 23	dim.SpatialImage
Centroids, 9, 25, 26, 30, 33, 53, 55, 56	(SpatialImage-methods), 107
Centroids-class, 24	dimensional reduction, 43, 84, 124, 128,
Centroids-methods, 25	129, 131
CheckGC, 27	dimensional reductions, 89
CheckLayersName, 27	dimnames.Assay, 14, 15, 22, 30, 42, 45, 46,
Command, 28, 81	47, 82, 111, 114, 123, 126, 132

dimnames.Assay5, <i>16</i> , 22, 42, 45, 46, 47, 83,	Graph, 9, 18, 105
112, 115, 124, 128, 133	Graph (Graph-class), 59
dimnames.DimReduc (dim.DimReduc), 43	Graph-class, 59
dimnames. Seurat, 22, 44–46, 46, 85, 89, 105,	Graphs (Assays), 18
117, 129, 131, 132, 135	grob, <i>109</i>
dimnames <assay (dimnames.assay),="" 45<="" td=""><td>hand Annua (FF Annua) 106</td></assay>	hand Annua (FF Annua) 106
dimnames <assay5 (dimnames.assay5),="" 46<="" td=""><td>head. Assay ([[. Assay), 126</td></assay5>	head. Assay ([[. Assay), 126
dimnames <seurat (dimnames.seurat),="" 46<="" td=""><td>head. Assay5 ([[.Assay5), 127</td></seurat>	head. Assay5 ([[.Assay5), 127
DimReduc, 18, 31, 39, 43, 84, 93, 104, 115,	head. Seurat ([[.Seurat], 129
116, 124, 128	here, 129, 131, 132
DimReduc (DimReduc-class), 47	HVFInfo, 59
DimReduc-class, 47	Idonto 62 104
DimReduc-validity, 48	Idents, 62, 104
Distances, 49	Idents<- (Idents), 62
drop, 124, 126–129	Images, 65
droplevels.LogMap, 11, 50, 67, 75, 80, 81	images, 89
droplevels. Seurat (Idents), 62	Index, 66
, , , , , ,	Index<- (Index), 66
Embeddings, 51, <i>128</i>	Indices, 66
	Inf, 9, 30, 33
Features (Cells), 22	intersect.LogMap, 11, 50, 67, 75, 80, 81
Features.FOV (FOV-methods), 54	is.finite.Centroids
Features. Molecules (Molecules-methods),	(Centroids-methods), 25
88	is.infinite.Centroids
FetchData, 51, 63, 122	(Centroids-methods), 25
FetchData.FOV (FOV-methods), 54	IsGlobal, 68, <i>110</i>
FilterObjects, 53	IsGlobal.SpatialImage
FOV, 33, 40, 54, 55	(SpatialImage-methods), 107
FOV (FOV-class), 53	IsMatrixEmpty, 68
FOV-class, 53	IsNamedList, 69
FOV-methods, 54	
FOV-validity, 57	JackStrawData, 70, 72
FOVs, 89	JackStrawData(JackStrawData-class), 70
1043, 67	JackStrawData-class, 70
GetAssayData (AssayData), 16	JackStrawData-methods, 70
GetImage, 57, 109, 110	JoinLayers, 71
GetImage.SpatialImage	JS, 72
(SpatialImage-methods), 107	JS<- (JS), 72
GetTissueCoordinates, 37, 58, 110	
GetTissueCoordinates.Centroids	Key, 73, 109, 110
	key, <i>112</i>
(Centroids-methods), 25	Key.SpatialImage
GetTissueCoordinates.FOV (FOV-methods),	(SpatialImage-methods), 107
54	Key<- (Key), 73
GetTissueCoordinates.Molecules	Key <spatialimage< td=""></spatialimage<>
(Molecules-methods), 88	(SpatialImage-methods), 107
GetTissueCoordinates.Segmentation	Keys (Key), 73
(Segmentation-methods), 101	Keys.FOV (FOV-methods), 54
GetTissueCoordinates.SpatialImage	
(SpatialImage-methods), 107	labels.LogMap, 11, 50, 67, 74, 80, 81

LayerData, 18, 75, 122-124	Neighbor, 11, 12, 18, 90
LayerData<- (LayerData), 75	Neighbor (Neighbor-class), 90
Layers (LayerData), 75	Neighbor-class, 90
layout, 58	Neighbor-methods, 90
length.Centroids (Centroids-methods), 25	Neighbors (Assays), 18
length.DimReduc(dim.DimReduc), 43	
length.FOV (FOV-methods), 54	Overlay, 91
lengths.Centroids(Centroids-methods),	Overlay, Centroids, SpatialPolygons-method (Overlay), 91
lengths.Segmentation	Overlay, FOV, FOV-method (Overlay), 91
(Segmentation-methods), 101	Overlay, FOV, Spatial-method (Overlay), 91
levels.Seurat (Idents), 62	Overlay, FOV, SpatialPolygons-method
levels <seurat (idents),="" 62<="" td=""><td>(0verlay), 91</td></seurat>	(0verlay), 91
Loadings, 77, 124	Overlay, Molecules, Spatial Polygons-method
Loadings<- (Loadings), 77	(Overlay), 91
LoadSeuratRds (SaveSeuratRds), 99	Overlay, Segmentation, Spatial Polygons-method
logical map, 11, 50, 67, 74	(Overlay), 91
logical mapping, 15	
LogMap, 11, 50, 67, 75, 78, 81, 106, 114	PackageCheck, 92
LogMap-class (LogMap), 78	pbmc_small, 92
LogMap-validity, 80	plan, <i>103</i>
LogSeuratCommand, 7, 10, 81, 106, 125, 135	plotly::layout, 58
	<pre>print (print.DimReduc), 93</pre>
Matrix, 9	print.DimReduc, 31, 43, 48, 49, 84, 93, 116,
matrix, 9, 36	124, 128
mean, <i>64</i>	Project, 36, 85, 94
merge (merge.Seurat), 84	Project<- (Project), 94
merge.Assay, 14, 15, 30, 42, 45, 82, 111, 114,	
123, 126, 132	Radius, 95, <i>110</i>
merge.Assay5, 16, 42, 46, 83, 112, 115, 124,	Radius.Centroids (Centroids-methods), 25
128, 133	Radius.SpatialImage
merge.DimReduc, 31, 43, 48, 49, 84, 94, 116,	(SpatialImage-methods), 107
124, 128	RandomName, 95
merge. Seurat, 44, 47, 84, 89, 105, 117, 129,	readRDS, 100
131, 132, 135	Reductions (Assays), 18
MergeSeurat (merge.Seurat), 84	regular expression, 77
Misc, 86	remove-object ([[<-,Seurat,NULL), 131
Misc<- (Misc), 86	remove-objects ([[<-,Seurat,NULL), 131
Molecules, 33, 34, 53, 55, 88	RenameAssays, 96
Molecules (Boundaries), 20	RenameCells, 97, 110
Molecules-class, 87	RenameCells.Centroids
Molecules-methods, 88	(Centroids-methods), 25
	RenameCells.FOV (FOV-methods), 54
name of a subobject, 129	RenameCells.Segmentation
names.DimReduc (dim.DimReduc), 43	(Segmentation-methods), 101
names.FOV (FOV-methods), 54	RenameCells.SpatialImage
names. Seurat, 44, 47, 85, 89, 105, 117, 129,	(SpatialImage-methods), 107
131, 132, 135	RenameIdent (Idents), 62
nearest-neighbor graphs, 89, 129	RenameIdents (Idents), 62

ReorderIdent (Idents), 62	SpatialPoints, 87
rlang::%, 104	split. Assay, 14, 15, 30, 42, 45, 82, 110, 114,
rlang::check_installed, 92	123, 126, 132
rle, 26, 103	
RowMergeSparseMatrices, 98	split.Assay5, 16, 42, 46, 83, 111, 115, 124, 128, 133
Rowner gespar seriati ices, 98	
sample, 95	StashIdent (Idents), 62
save, 99	Stdev, 113
saveRDS, 100	stdout, 71, 90, 106, 109
SaveSeuratRds, 99	StitchMatrix, 113
Segmentation, 9, 33–35, 53, 55, 56, 101, 102	subset, 110
Segmentation-class, 101	subset (subset. Seurat), 116
Segmentation-methods, 101	subset. Assay, 14, 15, 30, 42, 45, 82, 111,
set-if-null, 103	114, <i>123</i> , <i>126</i> , <i>132</i>
SetAssayData (AssayData), 16	subset. Assay5, 16, 42, 46, 83, 112, 115, 124,
SetDimReduction (CreateDimReducObject),	128, 133
31	subset. Centroids (Centroids-methods), 25
SetIdent (Idents), 62	subset.DimReduc, 31, 43, 48, 49, 84, 94, 115,
Seurat, 12, 19, 36, 39, 40, 44, 46, 47, 53, 85,	124, 128
89, 99, 116, 129–131, 134	subset.FOV (FOV-methods), 54
Seurat (Seurat-class), 104	subset. Molecules (Molecules-methods), 88
Seurat-class, 104	subset.Segmentation
Seurat-validity, 105	(Segmentation-methods), 101
	subset. Seurat, 44, 47, 85, 89, 105, 116, 129,
SeuratAccess (AddMetaData), 7	131, 132, 135
SeuratCommand, 7, 10, 81, 125, 135	subset.SpatialImage
SeuratCommand (SeuratCommand-class), 105	(SpatialImage-methods), 107
SeuratCommand-class, 105	SVFInfo (HVFInfo), 59
SeuratObject (SeuratObject-package), 6	toil Access (FF Access) 126
SeuratObject-package, 6	tail.Assay([[.Assay), 126
show, Centroids-method	tail.Assay5 ([[.Assay5), 127
(Centroids-methods), 25	tail.Seurat ([[.Seurat), 129
show, FOV-method (FOV-methods), 54	Theta, 117
show, JackStrawData-method	Theta.Centroids (Centroids-methods), 25
(JackStrawData-methods), 70	Tool, 105, 118
show, LogMap-method, 106	Tool<- (Tool), 118
show, Molecules-method	Tools (Tool), 118
(Molecules-methods), 88	UndataSauratObject 110
show, Neighbor-method	UpdateSeuratObject, 119 UpdateSlots, 120
(Neighbor-methods), 90	opuatestots, 120
show, Segmentation-method	v3, 89, 129, 131
(Segmentation-methods), 101	v5, 89, 129, 131
show, SpatialImage-method	valid0bject, 14-16, 48, 57, 80, 81, 105
(SpatialImage-methods), 107	VariableFeatures (HVFInfo), 59
Simplify, 106	VariableFeatures<- (HVFInfo), 59
SpatialImage, 53, 107	Version, 120
SpatialImage (SpatialImage-class), 107	voi 31011, 120
SpatialImage-class, 107	WhichCells, <i>116</i> , <i>117</i> , 121
SpatialImage-methods, 107	with_progress, 100, 103, 112
SnatiallyVariableFeatures (HVFInfo) 59	—, · · · · · · · · · · · · · · · · · · ·